



Instruction Manual

ZEISS Crossbeam 550L, Crossbeam 550

Focused Ion Beam - Scanning Electron Microscope (FIB-SEM)



ZEISS Crossbeam 550L, Crossbeam 550

Original Manual

Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena
Germany
info.microscopy.de@zeiss.com
www.zeiss.com/microscopy



Carl Zeiss Microscopy GmbH
Carl-Zeiss-Straße 22
73447 Oberkochen
Germany

Document Name: ZEISS Crossbeam 550L, Crossbeam 550 Instruction Manual
Order Number: 349500-8122-000
Revision: 3
Language: en-US
Effective from: 06/2020

© 2020 Without the prior written consent of ZEISS, this document or any part of it must neither be translated nor reproduced or transmitted in any form or by any means - including electronic or mechanic methods, by photocopying, recording or by any information or filing system. The right to make backup-copies for archiving purposes shall remain unaffected thereby. Any violations may be prosecuted as copyright infringements.

The use of general descriptive names, registered names, trademarks, etc. in this document does not imply that such names are exempt from the relevant intellectual property laws and regulations and therefore free for general use. This shall also apply if this is not specifically referred to. Software programs shall entirely remain the property of ZEISS. No program or subsequent upgrade thereof may be disclosed to any third party, copied or reproduced in any other form without the prior written consent of ZEISS, even if these copies or reproductions are destined for internal use at the customer's only, the only exception being one single back-up copy for archiving purposes.

ZEISS reserves the right to make modifications to this document without notice.

Table of Contents

| | |
|--|-----------|
| Glossary | 6 |
| 1 General Information | 9 |
| 1.1 Text Conventions and Link Types | 9 |
| 1.2 Explanation on Warnings and Additional Information | 10 |
| 1.3 Further Applicable Documents | 11 |
| 1.4 Contact..... | 12 |
| 2 Safety | 13 |
| 2.1 Intended Use | 13 |
| 2.2 Safety Procedures | 14 |
| 2.2.1 Safe Operation Conditions..... | 14 |
| 2.2.2 Requirements for Personnel..... | 14 |
| 2.2.3 Preventive Maintenance | 14 |
| 2.2.4 Safe Handling of Spare Parts | 15 |
| 2.3 Prevention of Hazards..... | 15 |
| 2.3.1 Biological Hazards | 15 |
| 2.3.2 Burn Hazards..... | 15 |
| 2.3.3 Chemical Hazards..... | 15 |
| 2.3.4 Electrical Hazards | 16 |
| 2.3.5 High Pressure Hazards..... | 16 |
| 2.3.6 Magnetic Field Hazards | 16 |
| 2.3.7 Mechanical Hazards | 16 |
| 2.3.8 Radiation Hazards | 17 |
| 2.3.9 Suffocation Hazards | 17 |
| 2.4 Safety Equipment..... | 18 |
| 2.4.1 Protective Cover Panels | 18 |
| 2.4.2 Main Switch | 19 |
| 2.4.3 Emergency Off (EMO) Option | 19 |
| 2.4.4 Locking Devices..... | 20 |
| 2.4.5 Lockout/Tagout Equipment | 21 |
| 2.4.6 Safety Labels | 22 |
| 3 Product and Functional Description | 27 |
| 3.1 System Overview | 27 |
| 3.2 Main Components | 29 |
| 3.2.1 Vacuum System..... | 29 |
| 3.2.2 Electron Optical Column Gemini II | 31 |
| 3.2.3 Ion-sculptor Focused Ion Beam (FIB) Column (Optional)..... | 35 |
| 3.2.4 Detectors | 37 |
| 3.2.5 Specimen Stage..... | 42 |
| 3.2.6 Dual Joystick | 43 |
| 3.3 Optional Components and Accessories | 44 |
| 3.3.1 Optional Detectors | 44 |
| 3.3.2 Control Panel | 53 |
| 4 Software Description | 55 |
| 4.1 SmartSEM..... | 55 |
| 4.1.1 SmartSEM User Interface..... | 55 |

| | | |
|----------|---|-----------|
| 4.1.2 | Graphical Control Elements | 56 |
| 4.1.3 | User Access Levels and User Privileges | 57 |
| 4.1.4 | SmartSEM Program Suite..... | 59 |
| 4.2 | SmartFIB..... | 61 |
| 4.2.1 | SmartFIB User Interface | 61 |
| 4.2.2 | Graphical Control Elements | 62 |
| 4.2.3 | SmartFIB Program Suite..... | 63 |
| 5 | Installation..... | 64 |
| 6 | Operation | 65 |
| 6.1 | Starting the System..... | 65 |
| 6.1.1 | Energizing the Microscope | 65 |
| 6.1.2 | Starting the Microscope | 66 |
| 6.1.3 | Starting SmartSEM | 66 |
| 6.2 | Obtaining a First Image..... | 68 |
| 6.2.1 | Preparing the Specimen Holder | 69 |
| 6.2.2 | Loading the Specimen Chamber | 70 |
| 6.2.3 | Locating the Specimen | 72 |
| 6.2.4 | Switching on the Gun..... | 73 |
| 6.2.5 | Switching on the EHT..... | 74 |
| 6.2.6 | Acquiring an Image | 74 |
| 6.2.7 | Optimizing the Image..... | 76 |
| 6.2.8 | Saving the Image..... | 78 |
| 6.3 | Modifying Gun Parameters | 79 |
| 6.3.1 | Selecting the Gun Mode..... | 79 |
| 6.3.2 | Setting the Probe Current..... | 80 |
| 6.3.3 | Measuring the Probe Current | 80 |
| 6.3.4 | Changing the Extractor Voltage..... | 81 |
| 6.4 | Working with Different Aperture Configurations and Beam Modes | 81 |
| 6.4.1 | Determining the Installed Aperture Configuration Gemini II Column | 81 |
| 6.4.2 | Selecting the Column Mode | 82 |
| 6.5 | Finding Appropriate Detector Settings | 82 |
| 6.5.1 | Selecting a Detector | 82 |
| 6.5.2 | Setting up the InLens SE Detector..... | 83 |
| 6.5.3 | Setting up the SE Detector | 83 |
| 6.5.4 | Setting up the EsB Detector..... | 84 |
| 6.5.5 | Setting up the SESI Detector..... | 85 |
| 6.5.6 | Setting up the aBSD/BSD Detector..... | 86 |
| 6.5.7 | Setting up the aSTEM/STEM Detector | 88 |
| 6.5.8 | Setting up the CL Detector | 91 |
| 6.6 | Shutting down the System..... | 91 |
| 6.6.1 | Finishing the Work Session | 91 |
| 6.6.2 | De-energizing the Microscope..... | 92 |
| 6.6.3 | Performing a Lockout/Tagout..... | 93 |
| 6.7 | Performing an Emergency Shutdown | 93 |
| 7 | Maintenance and Repair..... | 94 |
| 7.1 | Maintenance Work | 94 |
| 7.2 | Maintenance Intervals..... | 94 |
| 7.3 | Change of Consumables and Chemicals | 94 |

| | | |
|-----------|--|------------|
| 8 | Troubleshooting | 96 |
| 8.1 | Overview | 96 |
| 8.2 | Overall System | 99 |
| 8.2.1 | Checking the CAN Communication | 99 |
| 8.3 | Chamber | 99 |
| 8.3.1 | Initializing the Stage | 99 |
| 8.3.2 | Defining the Post Initialization Position of the Stage | 100 |
| 8.3.3 | Changing the Joystick TV Angle | 100 |
| 8.3.4 | Replacing the Chamber Door Seal | 101 |
| 8.3.5 | Resetting the Touch Alarm | 101 |
| 8.3.6 | Checking the Water Flow and Temperature | 102 |
| 8.4 | Column | 102 |
| 8.4.1 | Baking out the Gun Head | 102 |
| 8.4.2 | Calibrating the Probe Current | 103 |
| 8.5 | Power Circuit | 104 |
| 8.5.1 | Checking the Position of the Circuit Breakers | 104 |
| 8.6 | Detectors | 105 |
| 8.6.1 | Lubricating the Rod | 105 |
| 8.7 | Focused Ion Beam | 106 |
| 8.7.1 | Checking the Lifetime of the Ion Source | 106 |
| 8.7.2 | Regenerating the Ion Source by Heating | 106 |
| 8.7.3 | Manually Finding a New Extractor Value | 107 |
| 9 | Shutdown and Disposal | 108 |
| 9.1 | Putting the Microscope out of Operation | 108 |
| 9.2 | Transport and Storage | 108 |
| 9.3 | Disposal | 110 |
| 9.3.1 | Disposing of Solid Waste (Consumables) | 110 |
| 9.3.2 | Disposing of the Microscope | 110 |
| 10 | Technical Data and Conformity | 111 |
| 10.1 | Product Specification | 111 |
| 10.2 | Installation Requirements | 114 |
| 10.2.1 | Layout and Connections | 124 |
| 10.2.2 | System Layout | 127 |
| 10.3 | Declaration of Conformity | 129 |
| 11 | Parts and Tools | 130 |
| 11.1 | Consumables | 130 |
| 11.2 | Spare Parts | 130 |
| 11.3 | Tools and Accessories | 130 |
| | Index | 132 |

Glossary

aBSD

Annular Backscattered Electron Detector

AIC

Ampere Interrupting Capacity

Aperture

Mechanical limitation of an opening oriented perpendicular to the optical axis, which filters out electrons whose trajectories (tracks) do not run close to the optical axis.

aSTEM

Annular Scanning Transmission Electron Microscopy

Astigmatism

Lens aberration that distorts the shape of the electron beam, compensated by the stigmator.

Backscattered electrons

High-energy electrons that are liberated from the specimen surface when the specimen is hit by the primary electron beam.

Bakeout

Degassing of surfaces of a vacuum system by heating during the pumping process.

Beam booster

Anode and liner tube of the Gemini column are connected mechanically and electrically forming the beam booster. A booster voltage (UB, liner voltage) of +8 kV is applied to the beam booster, so that a high beam energy is maintained throughout the entire column. The beam booster technique has two main advantages: It minimizes beam widening, that may occur due to stochastic electron-electron interactions. Consequently there is almost no loss in beam brightness, even at low acceleration voltages. Secondly, the beam booster technique enhances protection against external stray fields.

BSD

Backscattered Electron Detector

BSE

Backscattered Electron

CAN

Controller Area Network

CCD

Charge-Coupled Device

CL

Cathodoluminescence

Condenser

Device that collects and focuses the electron beam onto the specimen.

D

Depth

DECT

Digital Enhanced Cordless Telecommunications

EBSD

Electron Backscatter Diffraction

EC

European Community

EDS

Energy Dispersive X-ray Spectroscopy

EDX

Energy Dispersive X-ray Spectroscopy

EHT

Extra High Tension

EIGA

European Industrial Gases Association

EM

Electron Microscope

EM server

A server that implements the internal communication between control software and microscope hardware.

| | |
|---|--|
| EMC | IGC |
| Electromagnetic Compatibility | Industrial Gases Council |
| EMO | IGP |
| Emergency Off | Ion Getter Pump |
| EsB | IR |
| Energy-selective Backscattered | Infrared |
| Eucentric | LOTO |
| Type of stage, the rotation axes of which intersect in the same point. The specimen surface is located in the eucentric point, where the tilt axis meets the beam axis. This guarantees that the focus is maintained when the specimen is tilted at a certain working distance. | Lockout/Tagout |
| Extractor | M |
| Positive electrode that attracts electrons from the filament. | M-axis |
| Faraday cup | MSDS |
| Small insulated metal container, equipped with an aperture where electrons can enter but not escape. Used to measure the specimen current in the microscope. | Material Safety Data Sheet |
| FESEM | PC |
| Field Emission Scanning Electron Microscope | Personal Computer |
| FIB | PE |
| Focused Ion Beam | Protective Earth (ground) |
| Focus wobble | PE |
| Function that sweeps the focus of the objective lens backwards and forward through the focus on the specimen plane. When the aperture is misaligned a lateral shift is observed. | Primary Electron |
| FTP | Penning gauge |
| File Transfer Protocol | Device for measuring high vacuum in the vacuum system. |
| GIS | Pre-vacuum pump |
| Gas Injection System | A pump for generating a pre-vacuum. |
| GUI | Primary electron beam |
| Graphical User Interface | Narrowly bundled beam of accelerated electrons that hit the specimen surface. |
| H | R |
| Height | R-axis (Rotation) |
| | RF |
| | Radio Frequency |
| | RMS |
| | Root Mean Square |
| | Schottky field emitter |
| | Type of electron source in which emission occurs at or near the work function barrier. |

Scintillator

Substance that absorbs electrons and in response, fluoresces photons while releasing the previously absorbed energy.

SE

Secondary Electron

Secondary electrons

Low-energy electrons that are emitted from the specimen surface when the specimen is hit by the primary electron beam. Secondary electrons are generated by inelastic scattering.

SEM

Scanning Electron Microscope

SEMI

Semiconductor Equipment and Materials International (SEMI) is a industry association comprising companies involved in the electronics design and manufacturing supply chain.

SESI

Secondary Electron Secondary Ion

SIMS

Secondary Ion Mass Spectrometry

STEM

Scanning Transmission Electron Microscope

Stigmator

Compensates astigmatism (lens aberration), so that the electron beam becomes rotationally symmetrical.

Suppressor

Electrode (anode) that suppresses unwanted thermionic emission from the shank of the Schottky field emitter.

T

T-axis (Tilt)

TEM

Transmission Electron Microscope

U

Voltage

User

Person examining a sample under the microscope.

W

Width

WD

Working Distance

WDS

Wavelength Dispersive X-ray Spectroscopy

WDX

Wavelength Dispersive X-ray Spectroscopy

WEEE

Waste Electrical and Electronic Equipment

X

X-axis

XeF₂ precursor

Xenon difluoride precursor

X-ray

Type of high energy electromagnetic radiation, that is generated during the operation of electron microscopes.

Y

Y-axis

Z

Z-axis

ZEISS Sales & Servicepartner

The sales and servicepartner is generally in the field for customer support in a regional area and / or a clearly defined customer group.

1 General Information

This manual is part of the SEM workstation, hereinafter referred to as the "microscope". Read the instructions carefully. Keep the manual nearby the microscope and hand it over to future owners of the microscope.

The manual is designed for operators who have been trained to operate the microscope by a ZEISS service representative. Basic operator training and safety instructions will be provided within the scope of initial start-up by ZEISS. Operators of the microscope must not deviate from the instructions provided in this manual.

This manual contains the following chapters:

| Chapter | Content |
|-------------------------------|--|
| Glossary | Lists important technical terms. |
| General Information | Explains the function and structure of this manual. |
| Safety | Summarizes important safety details. |
| Device Description | Describes the microscope and its main hardware components. |
| Software Description | Provides an overview of the user interface. |
| Installation | Refers to the ZEISS service representatives. |
| Operation | Contains information about starting the microscope and the software, obtaining a first image, adjusting important parameters, and powering down the microscope, also in emergency. |
| Maintenance and Repair | Informs you about preventive maintenance work and intervals and the change of consumables. |
| Troubleshooting | Describes common issues and how to resolve them. |
| Shutdown and Disposal | Summarizes notes on shutdown and disposal. |
| Technical Data and Conformity | Lists hardware specifications as well as the manufacturer's declaration that the equipment is in conformity with all applicable European Directives. |
| Parts and Tools | Lists consumables, spare parts, tools, and accessories. |
| Abbreviations | Lists abbreviations used in this manual. |
| Index | Lists keywords to help you find relevant information quickly. |

1.1 Text Conventions and Link Types

The following text conventions and link types are used:

| Text convention | Meaning |
|---|--|
| Click Start . Press the STANDBY button. Press [Enter] on the keyboard. | The names of controls and important information are shown in bold letters. |
| Press <Ctrl+Alt+Del> | Press several keys on the keyboard simultaneously. |

| Text convention | Meaning |
|---|---|
| Select Tools > Goto Control Panel > Air-lock . | Follow a path in the software. |
| <i>Text input</i> <i>Text output</i> | <ul style="list-style-type: none"> Text to be entered by the user Text displayed by the system |
| Programming and Macros | Anything typed in literally during programming, including, for example, macro codes, keywords, data types, method names, variables, class names, and interface names. |

Tab. 1: Text convention

| Link type | Meaning |
|---|---|
| See: <i>Text Conventions and Link Types</i> [▶ 9]. | Link to further information for this topic. |
| https://www.zeiss.com/corporate/int/home.html | Link to a website on the internet. |

Tab. 2: Link types

1.2 Explanation on Warnings and Additional Information

DANGER, WARNING, CAUTION, and NOTICE are standard signal words used to determine the levels of hazards and risks of personal injury and property damage. Not only the safety instructions and warnings in the **Safety** chapter are to be considered but also the safety instructions and warnings in other chapters. Failure to comply with these instructions and warnings can result in both personal injury and property damage and involve the loss of any claims for damages.

The following symbols and warnings indicating dangerous situations and hazards are used in this document.

DANGER

Type and source of danger

DANGER indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Type and source of danger

WARNING indicates a potentially hazardous situation which, if not avoided, may result in death or serious injury.

CAUTION

Type and source of danger

CAUTION indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

Type and source of danger

NOTICE indicates a potentially harmful situation which, if not avoided, may result in property damage.

Info

Provides additional information or explanations to help users better understand the contents of this Instruction Manual.

1.3 Further Applicable Documents

SmartSEM Software Manual

For detailed information on how to use the SmartSEM software for imaging and hardware control, refer to the SmartSEM software manual.

SmartFIB Software Manual

For detailed information on how to use the SmartFIB software for ion beam exposure, refer to the SmartFIB software manual.

Instruction Manual for Options

For detailed information regarding optional accessories, refer to the respective instruction manual in your document folder.

Product Specification and Installation Requirements

For details on technical data, refer to the documents Product Specification and Installation Requirements.

Material Safety Data Sheets

Material safety data sheets (MSDS) of chemicals used in combination with the microscope are contained in the document folder delivered with the microscope.

1.4 Contact

If you have any questions or problems, please contact your local ZEISS Sales & Servicepartner or one of the following addresses:

Headquarter

Phone: +49 1803 33 63 34

Fax: +49 3641 64 3439

Email: info.microscopy.de@zeiss.com

Service Germany

Phone: +49 7364 20 3800

Fax: +49 7364 20 3226

Email: service.microscopy.de@zeiss.com

Courses and training

Email: courses.microscopy.de@zeiss.com

ZEISS Sales & Servicepartner

You can find a ZEISS Sales & Servicepartner in your area under <https://www.zeiss.de/mikroskopie/website/forms/sales-and-service-contacts.html>.

2 Safety

This chapter contains general requirements for safe working practices. Any person using the Microscope System or commissioned with installation or maintenance must read and observe these general safety instructions. Knowledge on basic safety instructions and requirements is a precondition for safe and fault-free operation. Operational safety of the supplied Microscope System is only ensured if it is operated according to its intended use.

If any work is associated with residual risks, this is mentioned in the relevant parts of this document in a specific note. When components must be handled with special caution, they are marked with a warning label. These warnings must always be observed.

2.1 Intended Use

The microscope has been designed to generate an image or to analyze appropriate specimens, which is achieved by scanning a focused electron or ion beam across the specimen.

In addition to microscopic examination, the microscope also allows for the modification of appropriate specimens. For these purposes, the specimen is placed in the evacuated specimen chamber.

The microscope has been designed for the following applications:

- **Imaging**

Image generation and specimen analysis can be performed by means of a focused electron beam (SEM Imaging) or a focused ion beam (FIB Imaging) that is scanned across the specimen.

This application allows for the analysis of surface structures and near-surface structures of appropriate specimens.

- **Milling**

Milling comprises all applications that include material removal from the specimen surface by a focused ion beam.

In combination with a gas injection system option, the microscope has been designed for the following applications:

- **Gas-assisted deposition**

A process gas in combination with either a focused electron beam or a focused ion beam can be used to deposit material onto the specimen surface.

- **Gas-assisted etching**

This application enables an accelerated material removal by a focused electron beam or a focused ion beam in combination with a process gas.

Info

Not for therapeutic, treatment or medical diagnostic evidence.

Info

Not all products are available in every country. Contact your local ZEISS representative for more information.

SmartSEM Software

The SmartSEM software is intended for the operation of scanning electron microscopes (SEMs).

The SmartSEM software has to be run exclusively on a personal computer delivered by ZEISS. Any other applications are not allowed.

Regarding the operation of the microscope, the following regulations must be met:

- Only operate the microscope according to the operating conditions after correct installation by a ZEISS service representative.

- The microscope is only to be used by operators who have been trained by a ZEISS service representative. Basic operator training and safety instructions will be provided within the scope of the initial start-up by ZEISS. Make sure that everyone who operates the microscope only performs the tasks for which he/she is trained.
- Operators of the microscope must not deviate from the instructions provided in this manual.
- Only perform preventive maintenance tasks described in this manual. All tasks of maintenance, service, and repair not described in this manual have to be performed by an authorized ZEISS service representative.
- The microscope is to be used in a laboratory environment for commercial and scientific purposes only.

Using the microscope for any other purpose is not allowed and can be hazardous.

2.2 Safety Procedures

2.2.1 Safe Operation Conditions

If the product safety labels are covered or worn or if any of the safety devices are not in proper working condition, operation of the microscope can be hazardous.

- Periodically check the function of safety equipment.
- Periodically check that all protective cover panels are installed.
- Always follow the instructions given on the safety labels.
- Inspect and clean the product safety labels to maintain good legibility.

2.2.2 Requirements for Personnel

Deviating from the instructions given in this manual and on the safety labels can be hazardous or can lead to property damage.

- Do not operate the microscope until you have completely read and understood the entire user documentation delivered with the microscope.
- Observe all safety labels on the microscope and within this manual.
- Only operate the microscope according to the operating conditions after correct installation by a ZEISS service representative.
- Only ZEISS service representatives, who have specialized knowledge of radiation protection, are permitted to service the microscope.

Operator Training

Within the scope of initial start-up, the ZEISS service representative will perform a basic operator training. The basic operator training consists of fundamental operation procedures including safety instructions. Besides, an introduction to basic maintenance tasks will be given for an advanced operator, who has to be an electrically skilled person. The training performed must be documented appropriately.

Special application training is offered on request.

2.2.3 Preventive Maintenance

Deviating from the maintenance and repair tasks described in this manual can be hazardous or can lead to property damage.

- Only perform preventive maintenance and repair tasks described in this manual.
- All tasks of maintenance, service, and repair not described in this manual have to be performed by an authorized ZEISS service representative.

- To maintain best performance of the microscope, it is essential to perform preventive maintenance on a regular basis. Moreover, it is recommended that you conclude a service contract with your local ZEISS service representative or organization.

2.2.4 Safe Handling of Spare Parts

Using spare parts that are not provided by ZEISS can be hazardous or can lead to property damage:

- Only genuine parts supplied by ZEISS are to be used in servicing the microscope.
- Contact your ZEISS service representative for information regarding how to order spare parts.
- Unless authorized by ZEISS, all spare parts should be installed by a ZEISS service representative.

2.3 Prevention of Hazards

This section summarizes potential hazards and recommended safety precautions. Failure to follow the safety instructions and instructions may result in personal injury and property damage.

2.3.1 Biological Hazards

Danger to life: Biological substances may pose a threat to the health of humans and other living organisms.

Biological substances

- Keep a logbook of the biological substances loaded into the microscope and show it to the ZEISS service representatives before they perform any work on the microscope.

2.3.2 Burn Hazards

Risk of property damage: Hot surfaces during bakeout Parts of the enclosure in the upper range of the column may become hot during bakeout, particularly after a long bakeout cycle.

- Do not place any combustible objects on the grids of the electron optical column during bakeout.
- After the bakeout procedure, let surfaces cool down before working around the column.
- Only advanced operators are allowed to perform the bakeout procedure.

2.3.3 Chemical Hazards

Risk of injury: Aggressive or toxic chemicals Aggressive or toxic chemicals can cause chemical burns.

- When handling aggressive or toxic chemicals, wear suitable protective clothing, gloves, and eye/face protection.
- Do not eat, drink, or smoke while handling toxic chemicals.
- Refer to local safety regulations.
- Read and follow the instructions in the material safety data sheet of the chemical. The material safety data sheet can be obtained from the supplier of the chemicals.

Risk of injury: Reaction products Dangerous reaction products can be present in the specimen chamber during or after operation.

- Ensure that there is an appropriate exhaust gas line to remove the waste gas of the pre-vacuum pump and to transmit it to the outside.
- Wear lint-free gloves when touching the inner parts of the specimen chamber or the specimen.

Risk of injury: Irritant gases If a gas injection system (GIS) is used, irritant gases might be released from the precursors. Gases can cause irritation to eyes, skin, and respiratory system.

- Do not remove a reservoir from the gas injection system.
- Contact your ZEISS service representative to have an empty reservoir replaced.

- Never try to open a reservoir.
- For further information, refer to the GIS instruction manual.

Environmental risk: Disposal of aggressive or toxic chemicals

When disposing of aggressive or toxic chemicals, there is a threat of damage to the environment.

- When disposing of waste that has been generated during a service operation (e.g. used rotary pump oil), comply with all national and local safety and environmental protection regulations.

2.3.4 Electrical Hazards

Danger to life: Hazardous voltage inside the microscope

High voltages are present inside the microscope. Contact may cause burn or electrical shock.

- Ensure proper grounding. For more information, refer to the Installation Requirements document.
- Only authorized ZEISS service representatives are allowed to service the microscope.
- Do not try to service the microscope yourself.

Danger to life: High leakage currents

High leakage currents are present in the microscope. Contact may cause burn or electrical shock.

- Ensure proper grounding. For more information, refer to the Installation Requirements document.
- Do not operate the microscope without the separate ground connection.

2.3.5 High Pressure Hazards

Risk of injury or property damage: High pressure in gas cylinders

Gas cylinders containing for example nitrogen or compressed air have a high internal pressure of approximately 200 bar. If not properly handled, the contained gas can abruptly escape and cause the gas cylinder to move in an uncontrollable manner.

- Observe all safety labels on the gas cylinders and all safety instructions given by the gas cylinder manufacturer.
- Always operate gas cylinders in an upright position and secure them so they will not tip over.
- Before transporting gas cylinders, place protective caps on them.

2.3.6 Magnetic Field Hazards

Risk of injury: Malfunction of medical devices near ion getter pumps

Magnetic fields present at the ion getter pumps may disturb the function of medical devices. The magnetic fields are also present if the microscope is switched off.

If you wear medical implants that are susceptible to magnetic fields (e.g. cardiac pacemakers), do the following:

- Keep a distance of at least 1 m from the ion getter pumps.
- Follow the safety instructions provided by the pump manufacturer.

2.3.7 Mechanical Hazards

Risk of injury: Moving the specimen stage

Fingers can be trapped by the moving specimen stage.

- Always close the chamber door before moving the specimen stage.
- To remove parts fallen into or near to the stage use a tool (e.g. tweezers) instead of your fingers.

Risk of injury: Closing the chamber door

Fingers can be pinched when closing the chamber door.

- Use the recessed grip or door handle to close the chamber door.
- Ensure not to get your fingers caught in the chamber door gap.

- Risk of property damage: Short working distance** When opening the chamber door, the microscope or specimen can be damaged if the specimen stage is at a short working distance. If a BSD detector is inserted, it can be damaged as well.
- Always move the specimen stage to a long working distance before opening the chamber door.

2.3.8 Radiation Hazards

- Danger to life: X-rays** X-rays are generated inside the microscope during operation. This is unavoidable because electrons are accelerated by voltages up to 30 kV.
- Do not remove any parts around the column and chamber that are essential for radiation protection.
 - Use genuine ZEISS parts exclusively.
 - Ensure that all local safety and X-ray protection regulations are met.
 - Only authorized ZEISS service representatives are allowed to service the microscope.

The microscope is equipped with several radiation protection devices, which, under regular operating conditions, ensure that the microscope operates in accordance with the German X-ray protection regulation (RöV), the German radiation protection regulation (StrSchV) as well as with the EC Directive 2013/59/EURATOM.

In the EU, the operation of the microscope is permission-free as the following requirements are fulfilled:

- The acceleration voltage is limited to 30 kV.
- The local dose rate at a distance of 0.1 m from the accessible surface of the microscope does not exceed 1 µSv/h.
- A respective label is attached to the microscope.

Outside the EU, the user of the microscope has to comply with the local regulations of the country where the microscope is operated.

Radiation Protection

For questions regarding radiation protection, contact the ZEISS Radiation Safety Officer, Carl Zeiss AG, 73447 Oberkochen, Germany

phone: +49 (0) 7364 20 0

2.3.9 Suffocation Hazards

- Danger to life: Lack of oxygen** Gaseous dry nitrogen is used to vent the specimen chamber during specimen exchange and parts exchange. Inhaling nitrogen may cause unconsciousness.

- During specimen exchange, keep the chamber door open as short as possible.
- Do not inhale the air from within the specimen chamber.
- Ensure that the area around the microscope is sufficiently ventilated.
- If you begin to experience symptoms of asphyxia (for example: rapid breathing, loss of mental alertness and/or muscular coordination, depression of sensations, emotional instability, fatigue) leave the room immediately and inform the facility's safety officer.

Concerning the hazards of nitrogen installations and associated safety precautions refer to the current version of the guideline IGC Doc 44/18: Hazards of Oxygen-Deficient Atmospheres, published by EIGA (European Industrial Gases Association).

To download the document:

1. Go to EIGA homepage www.eiga.eu.
2. Select **Publications > EIGA Documents**.
3. From the list, select **Doc. 44/18**.
4. Click **Download**.

2.4 Safety Equipment

In order to prevent hazards to human health or property damage, the microscope is equipped with several safety features.

2.4.1 Protective Cover Panels

Due to hazardous voltages and X-rays inside the microscope, the microscope is equipped with protective cover panels.




Fig. 1: Protective cover panels

- 1** Electron optical column protective cover panels
- 2** Specimen chamber protective cover panels
- 3** Plinth protective cover panels

Operation of the microscope is only allowed with attached protective cover panels.

2.4.2 Main Switch

The Main Switch cuts off both phases of the mains power from the microscope.

The Main Switch guarantees an ampere interrupting capacity (AIC) of at least 10,000 A rms. In the OFF position, the Main Switch can be secured with a padlock against inadvertent activation, for example during repair and maintenance work.

The Main Switch is located at the rear side of the plinth.




Fig. 2: Main Switch

2.4.3 Emergency Off (EMO) Option

With the Emergency Off (EMO) option, the microscope is equipped with the following additional safety equipment.

Microscopes with SEMI certification must be equipped with the Emergency Off (EMO) option.

Emergency Off (EMO) Button

The EMO button is located on the plinth adjacent to the specimen chamber.

The EMO button must be pressed in an emergency to cut off mains power from the microscope and all devices connected to the AC Unit.

The EMO button must always be readily accessible and operable.

Up to three additional EMO buttons may be connected.



Fig. 3: EMO button

S2 (Start) Button

The green S2 (Start) button is attached below the MAIN switch at the rear side of the plinth.

The S2 (Start) button is required to confirm the setting of the MAIN switch. It cuts off all devices connected to the EMO box from the main power supply.

The S2 (Start) button must be pushed in the following cases:

- After the MAIN switch has been set to the ON position, or
- After the EMO button has been released.



Fig. 4: S2 (Start) button below the MAIN switch

2.4.4 Locking Devices

The microscope is equipped with several locking devices.

2.4.4.1 Chamber Door Locking Device

The chamber door locking device ensures that the door of the specimen chamber is closed properly.

It is located at the inner bottom front side of the specimen chamber.

When this locking device is activated (i.e. no electrical contact) `EHT Vac ready = no` is indicated in the SmartSEM user interface. EHT and SE detector voltages are blocked.

2.4.4.2 Vacuum Locking Device

The vacuum locking device ensures that gun vacuum and system vacuum are better than the required thresholds.

2.4.4.3 Interlock System of Optional Airlock

The interlock function ensures that the gate valve can only be operated if the chamber door and the airlock door are properly closed.




Fig. 5: Airlock Error LED

Additional blocking functions ensure that the specimen can only be transferred if the following conditions are fulfilled:

- The airlock rod is retracted.
- The specimen stage is in transfer position.
- The EHT (Extra high tension) is off.
- The column chamber valve is closed.

This is to prevent any risk of damaging the airlock rod or the gate valve.

2.4.5 Lockout/Tagout Equipment

The Lockout/Tagout (LOTO) equipment isolates the workstation from potential hazards associated with the unexpected release of hazardous energy.

Performing a Lockout/Tagout [▶ 93] can be necessary during

- maintenance tasks
- service work
- unscheduled repairs

At the site of installation, the electrical power connection and the fluid connections must be equipped with energy isolating devices.

These energy isolating devices (main shut-off devices)

- must be easily accessible
- must be mounted near the microscope in such a way that the person actuating or inspecting an energy isolating device is not be exposed to serious risks
- must close off the connections to the corresponding media when needed
- must be lockable in their off position in order to prevent accidental re-activation

The customer is responsible for

- having a working knowledge of the LOTO specification
- purchasing, distributing, and installing appropriate LOTO equipment
- auditing and enforcing compliance with LOTO procedures to all authorized personnel
- providing instructions on how to operate the energy isolating devices properly

Required Energy Isolating Devices

- Branch circuit molded case circuit breaker
- Main shut-off valves
 - Nitrogen supply
 - Compressed air supply
 - Cooling water supply


2.4.6 Safety Labels

Appropriate safety labels on the microscope warn operators of hazards. Each safety label is affixed close to the point where a particular hazard exists. Several labels also provide legal information.

Labels Attached to the Front Side of the Microscope




Labels Attached to the Rear Side of the Microscope



Labels Attached to the Gun Head underneath the Cladding










Labels Attached to the Rear Side of the Plinth






Labels Attached to the CEE Connector





The safety labels always need to be attached to the correct spots at the microscope. If a safety label is lost or unreadable, it needs to be reordered via the following reorder numbers:


| Position and Figure of the Safety Label | Description |
|--|--|
| <p>1 At the front of the chamber door.</p>  | <p>Risk of injury</p> <p>Fingers could be trapped. Always close the chamber door before you move the stage.</p> <p>Reorder no. 347800-0033-000-02en</p> |
| <p>2 At the front of the chamber door.</p>  | <p>Risk of damage</p> <p>The FESEM or specimen could be damaged if specimen stage is at a short working distance. Move the specimen stage to a long working distance before opening the chamber door.</p> <p>Reorder no. 347800-0033-000-04en</p> |
| <p>3 At the front of the plinth</p>  | <p>Avoid injury</p> <p>Make sure you have read and understood the instruction manual before operating this product.</p> <p>Reorder no. 347800-0033-000-01en</p> |
| <p>4 At the rear of the electron optical column</p>  | <p>Hazardous voltage inside</p> <p>Contact may cause burn or electric shock. Only authorized service Staff is allowed to service the equipment.</p> <p>Disconnect power before opening.</p> <p>Reorder no. 347800-0033-000-03en</p> |
| <p>5 At the ion getter pumps</p>  | <p>Magnetic field</p> <p>Can be harmful to pacemaker wearers. Pacemaker wearers stay back 30cm (12 in.).</p> <p>Reorder no. 360735-0000-034en</p> |
| <p>6 At the FIB column</p>  | <p>Hazardous voltage symbol</p> <p>Reorder no. 360400-0000-215</p> |
| <p>7 At the rear of the microscope</p>  | <p>Suffocation hazard</p> <p>The specimen chamber is ventilated with gaseous nitrogen.</p> <p>Ensure that the area around the electron microscope is sufficiently ventilated.</p> <p>Reorder no. 347800-0033-000-08en</p> |

| Position and Figure of the Safety Label | Description |
|---|---|
| <p>8 At the rear of the microscope</p>  | <p>High leakage current</p> <p>Ensure proper grounding.</p> <p>Do not operate the electron microscope without a separate ground connection.</p> <p>Reorder no. 347800-0033-000-07en</p> |
| <p>9 At the rear of the microscope</p>  | <p>Radiation hazard</p> <p>Do not remove any parts.</p> <p>Observe local safety and X-ray protection regulations.</p> <p>Service only through authorized personnel.</p> <p>Reorder no. 347800-0033-000-06en</p> |
|  | <p>For relevant options:</p> <p>Liquid nitrogen</p> <p>Contact with liquid or cold vapors can cause freeze burns or severe frostbite.</p> <p>Wear face, eye and skin protection.</p> <p>Reorder no. 347800-0033-000-09en</p> |

Underneath the cover panels of the microscope, there are some more safety labels, which are addressed to authorized ZEISS service representatives. These safety labels are described in the documents for ZEISS service representatives.

| Position and Figure of the Type Plate | Description |
|--|--------------------------|
| <p>At the rear of the microscope</p>  | <p>Type plate</p> |

| Position and Figure of the Legal Information | Description |
|--|--|
| <p>At the rear of the microscope</p>  | <p>Caution</p> <p>X-rays are produced in this Scanning Electron Microscope!</p> <p>The acceleration voltage is limited to 30 kV.</p> <p>Observe local safety and X-ray protection regulations.</p> <p>Order no. 346000-0088-000</p> |

| Position and Figure of the Type Plate | Description |
|--|---|
| <p>11 At the rear of the microscope</p>  | <p>The Type plate contains a unique serial number, the Sigma model designation with the according part number, the manufacturing year, and the mains voltage for which the Sigma is configured.</p> <p>Order no. 356100-9057-000 (Sigma 300) Order no. 352100-9903-000 (Sigma 500)</p> |

| Position and Figure of the Type Plate | Description |
|--|--|
| <p>13 At the rear of the microscope</p> | <p>The Type plate contains a unique serial number, the EVO model designation with the according part number, the manufacturing year, and the mains voltage for which the EVO is configured, either 230 VAC or 120 VAC.</p> <p>Order no. 354900-9040-000 (EVO10) Order no. 354900-9041-000 (EVO15) Order no. 354900-9042-000 (EVO25)</p> |

3 Product and Functional Description

3.1 System Overview

Main Components




Fig. 6: Main components of the Microscope System

- | | |
|--|--|
| <p>1 GEMINI II electron optical column Refer to <i>Electron Optical Column Gemini II</i> [▶ 31]</p> <p>3 Ion-sculptor focused ion beam (FIB) column Refer to <i>Ion-sculptor Focused Ion Beam (FIB) Column (Optional)</i> [▶ 35]</p> <p>5 Plinth with ON/STANDBY/OFF buttons</p> <p>7 Work desk</p> <p>9 Personal Computer (PC)</p> | <p>2 Ion getter pumps (IGPs) for FIB</p> <p>4 Specimen chamber with door handle</p> <p>6 Monitor</p> <p>8 Dual joystick Refer to <i>Dual Joystick</i> [▶ 43]</p> |
|--|--|

Optional Components and Accessories




Fig. 7: Optional components and accessories

- | | |
|--|---|
| <p>1 Multichannel Gas Injection System (GIS) for up to 5 precursor materials on a single flange Refer to instruction manual Multi GIS</p> | <p>2 Annular STEM detector (aSTEM) for transmission imaging and quality control. Two different models are available (aSTEM4, aSTEM1). Detection of up to 4 channels in parallel is possible. Refer to <i>aSTEM Detector</i> [▶ 50]</p> |
| <p>3 Local charge compensator for SEM imaging and analysis of non-conductive specimens Refer to instruction manual Charge Compensator</p> | <p>4 Electron flood gun for ion beam preparation of non-conductive specimens Refer to instruction manual Flood Gun</p> |
| <p>5 Single needle GIS for high angle specimen access Refer to instruction manual UniGIS</p> | <p>6 Micro-manipulator for specimen modification and probing</p> |
| <p>7 BSD detector for high efficiency BSE detection and angle selective material characterization. Four different models are available (aBSD4, aBSD1, BSD4, BSD1). Detection of up to 4 channels in parallel is possible. Refer to <i>BSD Detector</i> [▶ 46] and <i>aBSD Detector</i> [▶ 47]</p> | <p>8 EsB detector for finest z resolution without topographic artifacts and unique material contrast Refer to <i>EsB Detector</i> [▶ 45]</p> |
| <p>9 Airlock with optional navigation camera (80 mm or 200 mm wide) for fast and efficient specimen transfer and fast pumping times Refer to the respective instruction manual 80 mm Airlock or 200 mm Airlock</p> | |

Further options:

- Control panel
Refer to *Control Panel* [▶ 53]
- Plasma Cleaner
Refer to instruction manual Plasma Cleaner
- Energy-dispersive X-ray spectrometer (EDS) for chemical analysis
- EBSD for crystallographic mapping
- SESI detector for Secondary Electron and Secondary Ion imaging
Refer to *SESI Detector* [▶ 44]
- CL detector for the analysis of cathodoluminescent materials
Refer to *CL Detector* [▶ 52]
- Nanopatterning and Visualization Engine (NPVE) for advanced patterning and lithography tasks
Information will be provided by the manufacturer
- Electrostatic Beam Blanker for SEM
Refer to instruction manual Beam Blanker
- Rapid Laser Ablation Upgrade
Refer to instruction manual Rapid Laser Ablation Upgrade
- Further options on request

3.2 Main Components

3.2.1 Vacuum System

Purpose For operation of the microscope, the gun head, the column, and the specimen chamber have to be evacuated. The vacuum is essential to operate the gun and to prevent collisions of electrons with gas molecules.




Fig. 8: Schematics of the vacuum system of a microscope with Gemini II column

- 1** Gun with filament
- 2** Gun head
- 3** Single hole aperture
- 4** Column chamber valve
- 5** Turbo pump
- 6** Ion getter pumps (IGP)
- 7** Specimen chamber
- 8** Penning gauge
- 9** Vent valve
- 10** Pre-vacuum pump

Operating Principle

The pre-vacuum pump **10** and the turbo pump **5** evacuate the specimen chamber **7**. The vacuum in the specimen chamber is measured by a Penning gauge **8**. The detected vacuum values are displayed as System vacuum in the SmartSEM user interface. As long as the detected pressure in the specimen chamber is not ready for operation, the column chamber valve **4** is closed in order to separate the specimen chamber from the column.

Gun Vacuum

In the gun head, an ultra high vacuum is maintained by ion getter pumps **6**. The vacuum in the gun head is displayed as Gun vacuum in the SmartSEM software. It should be below 1×10^{-8} mbar.

Venting

The specimen is located in the evacuated specimen chamber. To open the specimen chamber for specimen exchange, you have to break the vacuum in a controlled manner. This is done by the Vent command via the SmartSEM user interface or by pressing the **Exchange** push button on the optional control panel.

When the Vent command is received, the column chamber valve closes and gaseous nitrogen flows into the specimen chamber via the vent valve **9**. As soon as the pressure equilibrium is obtained, the chamber door can be opened to change the specimen.

Evacuating

In order to continue operation, the Pump command makes the pre-vacuum pump and the turbo pump evacuate the specimen chamber.

As soon as the vacuum in the specimen chamber is ready for operation, the column chamber valve opens and the EHT Vac ready message is displayed in the SmartSEM user interface. Gun and EHT can be switched on.

Quiet Mode

The automatically controlled **Quiet Mode** is optionally available. This option allows switching off the pre-vacuum pump after specimen exchange when the vacuum threshold is achieved.

3.2.2 Electron Optical Column | Gemini II

Purpose The GEMINI II column is the part of the microscope, where electrons are emitted, accelerated, deflected, focused, and scanned. Main characteristics of the GEMINI optics are the beam booster and an objective lens that consists of a combined electrostatic/electromagnetic lens doublet.




Fig. 9: Schematics of the electron optics

- 1** Gun
- 2** Extractor
Positive electrode that extracts electrons from the filament.
- 3** Anode aperture
- 4** Upper condenser
- 5** Single hole aperture
- 6** Lower condenser
- 7** EsB detector
- 8** InLens SE detector
- 9** Objective lens
Focuses the electron beam on to the specimen surface.
- 10** Scanning coils
Deflect the beam across the specimen surface in what is usually referred to as a raster scan.
- 11** Specimen
- U_{Sup}** Suppressor voltage
- U_{Ext}** Extractor voltage
- U_{EHT}** Acceleration voltage
- U_B** Liner tube voltage

Operating Principle

A Schottky field emitter serves as gun **1**. The filament is heated by applying the filament current. Electrons are emitted from the heated filament while an electrical field, called extractor voltage (U_{Ext}), is applied. To suppress unwanted thermionic emission from the shank of the Schottky field emitter, a suppressor voltage (U_{sup}) is applied as well.

EHT

The emitted electrons are accelerated by the acceleration voltage (U_{EHT}). The beam booster (U_{B} , booster voltage), which is always at a high potential if the acceleration voltage is 20 kV or less, is integrated directly after the anode. This guarantees that the energy of the electrons in the entire beam path is always much higher than the set acceleration voltage. This considerably reduces the sensitivity of the electron beam to magnetic stray fields and minimizes the beam broadening.

Apertures

The electron beam passes through the anode aperture **3** first, afterwards through a single hole aperture **5**.

The anode aperture defines the maximum possible probe current.

For the Gemini II column, two different column configurations are available:

- 40 nA high resolution configuration

| Anode aperture diameter | Probe current | Typical application |
|-------------------------|----------------|---------------------|
| 55 μm^* | 10 pA to 40 nA | High resolution |

- 100 nA high current configuration

| Anode aperture diameter | Probe current | Typical application |
|-------------------------|-----------------|--|
| 90 μm^* | 10 pA to 100 nA | Combined high resolution and analytical investigations |

Condenser

A double condenser system allows the continuous regulation of the probe current. The upper condenser **4** is used to continuously adjust the probe current. The lower condenser **6** is used for aperture matching of the objective lens in order to guarantee optimum resolution at each probe current and EHT setting.

Stigmator

The stigmator is located inside the condenser and compensates for astigmatism so that the electron beam becomes rotationally symmetrical.

Deflection System

The electron beam is focused by the objective lens **9** onto the specimen **11** while being deflected in a point-by-point scan over the specimen surface by the scanning coils.

Before the electron beam exits the objective lens, the electrostatic lens creates an opposing field which reduces the potential of the electrons by +8 kV. The energy of the electrons reaching the specimen surface therefore corresponds to the set acceleration voltage (EHT).

Signal Detection


When the primary electron beam hits the specimen, certain interaction products are released, which can be recorded by specific detectors, e.g. the InLens SE detector **8**. For more information see *Principle of Signal Detection* [[▶ 37](#)].

3.2.2.1 Beam Modes

With the GEMINI II column, different beam modes are available:

- High Resolution
- High Resolution Low Current
- Normal
- Minimum Probe Current
- Maximum Probe Current
- Fisheye
- Widefield


High Resolution

| Characteristics | Beam Path |
|---|--|
| <p>The upper condenser 1 has a low excitation which is varied in order to adjust the probe current in a limited range.</p> <ul style="list-style-type: none"> ▪ Electron interaction effects are minimized ▪ High spatial resolution is guaranteed ▪ Probe current range is limited |  |

High Resolution Low Current

The characteristics and the beam paths of the High Resolution Low Current mode and High Resolution mode are the same, except the different operation of the gun.

Normal

| Characteristics | Beam Path |
|--|---|
| <p>The upper condenser 1 focuses the beam and creates a crossover above the beam-limiting aperture 2. The upper condenser is adjusted to modify the probe current.</p> <p>The lower condenser 3 is used for aperture matching of the objective lens 4 for optimum spatial resolution.</p> <ul style="list-style-type: none"> ▪ Whole probe current range is accessible ▪ High spatial resolution is usually guaranteed <p>Instead of aperture matching it is also possible to adjust the amount of depth of field desirable. The lower condenser 3 is then operated at increasing excitations.</p> <ul style="list-style-type: none"> ▪ Provides a higher depth of field ▪ Useful to investigate specimens with a high aspect ratio or to navigate on a tilted specimen |  |

Minimum Probe Current

The characteristics and the beam paths of Minimum Probe Current mode and Normal mode are the same, except the different operation of the gun.

Maximum Probe Current

The characteristics and the beam paths of Maximum Probe Current mode and Normal mode are the same, except the different operation of the gun.

Fisheye

- Provides a huge field of view, allows you to image the whole specimen holder with the electron beam.
- Useful for an overview of large specimens and for specimen navigation.

Widefield

- Provides a huge field of view, allows you to image the whole specimen holder with the electron beam.
- Useful for an overview of large specimens and for specimen navigation.
- Can be used with EHTs between 100 V and 2000 V.
- Requires the SmartSEM software license Field Mode.

3.2.3 Ion-sculptor Focused Ion Beam (FIB) Column (Optional)

Purpose The Ion-sculptor Focused Ion Beam (FIB) column is the part of the microscope, where ions are emitted, accelerated, focused, and deflected.




Fig. 10: Schematics of the ion optics

- 1** Ion source
Liquid metal ion source of gallium (Ga^+).
- 2** Variable apertures
- 3** Ion beam
- 4** Objective lens
- 5** Specimen

Operating Principle Gallium ions (Ga^+) are extracted from a liquid metal ion source **1**. The ions are accelerated by the acceleration voltage to an energy of maximum 30 keV. The ion emission is adjusted by the extractor and the suppressor voltage.

Gallium is used up during operation. Therefore, the gallium emitter cartridge is a consumable.

The gallium emitter has to be regenerated by heating from time to time; the heating procedure removes the gallium oxide, which has been created during operation.

Condenser

The electrostatic condenser collimates and focuses the ion beam.

Probe Currents

After passing the condenser, the beam current is defined by a set of software-controlled mechanical apertures. By using the different apertures in combination with the different condenser settings, the probe current can be continuously adjusted in the range between 1 pA and 100 nA.

Objective Lens

The objective lens is designed as an electrostatic Einzel-lens system. It focuses the beam onto the specimen surface.

3.2.3.1 Imaging Modes

Purpose With the FIB option, different imaging modes are available to operate the microscope.




Fig. 11: Imaging modes available with the FIB option

| Imaging Mode | FIB Mode | Characteristics | Typical Application |
|------------------------------|------------|--|--|
| 1 SEM imaging | SEM | <ul style="list-style-type: none"> ▪ Electron beam is active. ▪ Ion beam is blanked. ▪ The SE signal is synchronized to the SEM scan. | High resolution imaging. |
| 2 FIB imaging | FIB | <ul style="list-style-type: none"> ▪ Electron beam is blanked. ▪ Ion beam is active. ▪ The SE signal is synchronized to the FIB scan. | <ul style="list-style-type: none"> ▪ Channeling contrast imaging. ▪ Voltage contrast imaging. ▪ Defining milling patterns on the specimen surface. ▪ Grain analysis. |
| 3 Crossbeam operation | SEM + FIB | Image is composed of SEM and FIB components. | Setting the coincidence point. |
| | Mill | <ul style="list-style-type: none"> ▪ No image. ▪ Mills with the milling parameters set (milling current). | <ul style="list-style-type: none"> ▪ Ion beam milling or ion beam deposition. |
| | Mill + SEM | Mills and generates a SEM image. | SEM real-time imaging of the ion beam milling or deposition. |

3.2.4 Detectors

3.2.4.1 Principle of Signal Detection

The interaction products most frequently used for the generation of images in scanning electron microscopy are secondary electrons (SEs) and backscattered electrons (BSEs).

Primary Electrons (PEs)

Primary electrons are electrons forming the scanning beam before hitting the specimen.

Secondary Electrons (SEs)

Secondary electrons are emitted from the topmost layer of the specimen.

- **SE1 Electrons**

Electrons emitted at the point of impact between the beam and the specimen are known as SE1 type electrons. The amount of electrons emitted at the point of impact is related to the shape of the specimen.

Secondary electron detectors, such as the InLens SE detector, collect SE1 type electrons from the surface layer of the specimen and are thus ideal for displaying surface structures.

- **SE2 Electrons**

The emergence of backscattered electrons from the specimen excites further emission of secondary electrons. These are known as SE2 type electrons.

Detectors that collect SE2 type electrons are especially suitable where the working distance is large. Surface detail as the effect of "lateral illumination" emphasizes the topography of the specimen.

Backscattered Electrons (BSEs)

Backscattered electrons (BSEs) emerge from below the surface of the specimen (up to an order of μm). The number of electrons emitted at the point of impact is highly dependent on the mean atomic number of the specimen material. This means that a BSE image provides depth information and atomic number contrast.

BSE detectors are used to display the materials contrast because the backscatter coefficient is dependent on the mean atomic number of the material under investigation.

Transmitted Electrons

This comprises primary electrons that are transmitted through an ultrathin specimen and weakly scattered primary electrons with a small range of angles. Depending on the material, primary electrons are scattered under different angles and can be detected by a STEM detector placed below the specimen. Unscattered electrons are detected in the center of the STEM detector and give a bright field image. Electrons scattered under higher angles are detected by outer areas of the STEM detector and produce dark field images.

Cathodoluminescence (CL)

Electrons impacting on luminescent materials cause the emission of photons which may have wavelengths in the visible spectrum and can be imaged by specialized detectors.

3.2.4.2 Detectors Overview

The beam scans the specimen and initiates particles to be emitted. A detector collects the emission and produces an electric signal with an amplitude proportional to the number of particles at any given time.

| Standard Detectors | Detected Signals | Typical Application |
|--|-----------------------|---|
| <i>InLens SE Detector</i> [▶ 39] | SE1 | Surface structure |
| <i>SE Detector</i> [▶ 40] | SE2 | Topography |
| Optional Detectors | Detected Signals | Typical Application |
| <i>SESI Detector</i> [▶ 44] | SE2 Ions | Channeling contrast (crystal orientation) Compositional contrast |
| <i>EsB Detector</i> [▶ 45] | SE1 BSE | Surface structure Material contrast |
| <i>BSD Detector</i> [▶ 46] | BSE | Compositional contrast |
| <i>aBSD Detector</i> [▶ 47] | BSE | Angular resolved BSE imaging Compositional contrast |
| <i>aSTEM Detector</i> [▶ 50] | Transmitted electrons | Diffraction contrast Compositional contrast |
| <i>CL Detector</i> [▶ 52] | Light photons | Mineralogy |

3.2.4.3 InLens SE Detector

Purpose The InLens SE detector is a high-efficiency detector for high resolution SE imaging and detects secondary electrons directly in the beam path. The detection efficiency of this detector results from its geometric position in the beam path and from the combination with the electrostatic/ electromagnetic lens.

Position The annular shaped in-column detector is located above the objective lens.




Fig. 12: Schematics of the InLens SE detector

- 1 InLens SE detector
- 2 Beam path
- 3 Objective lens
- 4 Specimen

Operating Principle The primary electrons are accelerated by the acceleration voltage at the anode and, up to an acceleration voltage of 20 kV, by an additional beam booster voltage of 8 kV at the liner tube. To ensure that the electrons reach the specimen surface **4** with the energy set as acceleration voltage, an opposing electrostatic field of the same magnitude as the beam booster voltage (8 kV) is generated at the end of the objective lens by the electrostatic lens. This electrostatic field acts as acceleration field for the SEs generated on the specimen surface.

At the InLens SE detector, the electrons hit a scintillator. This generates a flash of light that is led out of the beam path and onto a photomultiplier by means of a lightguide. The photomultiplier converts the light information into an electronic signal, which can be displayed on the monitor.

The efficiency of the InLens SE detector is mainly determined by the electric field of the electrostatic lens, which decreases exponentially with the distance.

Thus, the working distance (WD) is one of the most important factors affecting the signal-to-noise ratio of the InLens SE detector.

As the tilt angle of the specimen surface affects the emission angle of the electrons, you should avoid strong specimen tilting.

Info

The InLens SE detector can be used up to an acceleration voltage of 20 kV. At higher acceleration voltages, the beam booster and thus the field of the electrostatic lens are switched off. Without the field of the electrostatic lens, which attracts the secondary electrons, the efficiency of the InLens SE detector is reduced.

3.2.4.4 SE Detector

Purpose The SE detector is an Everhart-Thornley type detector. It detects SEs as well as BSEs.

Position The SE detector is attached to the options ports plate.




Fig. 13: Schematics of the SE detector

- 1** Preamplifier
- 2** Photomultiplier
- 3** Light guide
- 4** Scintillator
- 5** Collector grid
- 6** Specimen

Operating Principle Electrons moving to the detector are attracted/repelled by the collector grid and directed to the scintillator. The electrons gain energy from the scintillator and thus are able to interact with a phosphor layer, which generates photons (light). The light travels up a light pipe to a photomultiplier. The photomultiplier multiplies the flashes of light and outputs a signal that can be used for imaging.

The collector voltage can be varied in the range between -250 V and $+400\text{ V}$.

A positive collector voltage generates an electrical field in front of the detector, thus directing the low energy SEs towards the scintillator.

Info

For all standard applications, the collector bias should be set to $+300\text{ V}$.

3.2.4.5 Chamber CCD Camera

Purpose The microscope contains a CCD camera (charge-coupled device camera) inside the specimen chamber. It is referred to as the chamber CCD camera or chamberscope. It allows you to monitor the position of the specimen stage and particularly the distance between the objective lens and the specimen holder.

The chamber CCD camera is located at the backside of the specimen chamber. A second CCD camera can be installed on the left-hand side of the specimen chamber.




Fig. 14: Sample image from chamber CCD camera.

- 1** Objective lens of the SEM column
- 2** Objective lens of the FIB column
- 3** Specimen holder

NOTICE

Risk of property damage: Risk of collision

Use the chamber CCD camera to monitor the position of the specimen holder during stage movements. Pay particular attention to the distance between the objective lens and the top of the specimen. This applies to vertical movements, but also to horizontal movements, because a thick specimen may collide with the objective lens from the side.

The chamber CCD camera has two illumination modes. The chamber can be illuminated either with white light or with infrared light. Infrared light gives a grayscale image, whereas white light gives a color image. In standard settings the mode is automatically selected, depending on the imaging mode and the selected detectors. White light limits the performance of most detectors. Therefore infrared illumination is a fallback if white light cannot be used. The performance of diode detectors is negatively affected by infrared light. If a diode detector is selected, then by default the chamber CCD camera is disabled. The automatic selection of the illumination mode can be manually overwritten by the user.

The image of the chamber CCD camera can be adjusted by the user:

- Brightness, contrast, and CCD illumination can be optimized to improve the image.
- Certain display details can be chosen by zooming in and shifting. This can be achieved by turning the Magnification and Shift X or Shift Y knobs on the control panel, respectively.

Info

It is highly recommended that the CCD illumination control is set to **Auto Detect**. This avoids any problems by the user forgetting to switch the illumination back for a different detector.



Fig. 15: Chamber CCD camera disabled as indicated by a pause sign (e.g. if a diode detector is selected).

3.2.5 Specimen Stage

Purpose The motorized eucentric specimen stage is used to navigate the specimen inside the specimen chamber.

Position The specimen stage is mounted on the chamber door. If the chamber door is closed, the specimen stage is inside of the specimen chamber.




Fig. 16: Specimen stage with dovetail fitting for precise fitting

- 1** Specimen holder
- 2** Specimen stage

Operating Principle When moving a tilted stage in X/Y direction, no displacement (moving out of focus) occurs. When a pre-tilted specimen holder is used, X/Y movements require an adjustment of the Z position to stay in focus.

The stage can be operated using the *Dual Joystick* [▶ 43] or using the SmartSEM software.

| Axis | Description | Movement |
|------|-------------|---|
| X | X-axis | Movement towards or away from the chamber door (horizontal movement in the image) |
| Y | Y-axis | Movement |
| Z | Height | Vertical movement (movement towards or away from the focal plane of the image) |
| M | Height | Movement of the specimen surface into the tilt axis at various working distances |
| R | Rotation | Stage rotation parallel to the X-Y plane |
| T | Tilt | Stage tilt about an axis parallel to the X axis |

The stage is eucentric, which means that the tilt axis intersects the beam axis. The specimen surface is located in the eucentric point, where the tilt axis meets the beam axis. This guarantees that the focus and the point of interest are maintained when the specimen is tilted at a certain working distance.

3.2.6 Dual Joystick

Purpose The dual joystick is used for stage control and specimen navigation.

Position The dual joystick is placed on the microscope desk.



Fig. 17: Dual joystick

1 T/Z joystick

Controls the Z axis and the stage tilt (T).

2 Break push button

Stops the stage in an emergency.

3 M push buttons

Controls a second Z-axis (M) on super-eucentric stages to set the eucentric point of the specimen tilt on these stages.

4 X/Y/R joystick

Controls the X- and Y-axes and the stage rotation by turning.

Operating Principle All axes are deflection-compensated. When the joystick is moved only slightly, the respective axis moves slowly. Larger movements of the joystick result in a faster movement of the stage.

The X-, Y-, and Z-axes are magnification-compensated. When working at a low magnification, the stage moves relatively fast. At higher magnifications the stage movement is slower. The stage is moving with its maximum speed when viewing the specimen with the CCD (Charge Coupled Device) camera.

The different axes can also be moved simultaneously.

3.3 Optional Components and Accessories

3.3.1 Optional Detectors

3.3.1.1 SESI Detector

Purpose The Secondary Electrons Secondary Ions (SESI) detector is suitable to detect secondary electrons as well as secondary ions. In Crossbeam systems, the optional SESI detector replaces the SE detector.




Fig. 18: Schematics of the SESI detector

- 1** FIB column
- 2** Objective lens
- 3** SESI detector
- 4** Collector grid
- 5** Specimen

Operating Principle Depending on the polarity of the collector voltage, either electrons or ions scattered from the specimen **5** are attracted by a collector grid **4** and accelerated to the converter. In the converter, both electrons and ions are converted into secondary electrons which are used to generate an image.

| Detector Mode | FIB Mode | Detected Signals | Typical Application |
|---|----------|---------------------|---|
| SE mode typical collector voltage +300 V | SEM, FIB | Secondary electrons | Topography |
| Ion mode typical collector voltage -4 kV | FIB | Secondary ions | Crystal orientation contrast, material contrast e.g. imaging of corrosion/oxidation caused by FIB processes in metals |

3.3.1.2 EsB Detector

Purpose The Energy-selective Backscattered (EsB) detector enables both SE and BSE detection. It can be used for high contrast topography or for compositional contrast examinations.

The SEs and BSEs generated at the impact point of the primary electron beam are intercepted by the low electrical field of the column. These electrons are accelerated by the field of the electrostatic lens.

Position The annular shaped in-column detector is located above the InLens SE detector.




Fig. 19: Schematics. SE detection via InLens SE detector (left) and BSE detection via EsB detector (right)

- 1 EsB detector
- 2 Filtering grid
- 3 InLens SE detector

Operating Principle

A small amount of SEs pass through the hole of the InLens SE detector **3** and would be observed by the EsB detector **1**. To prevent detection of these SEs, a filtering grid **2** is installed in front of the EsB detector. By switching on the filtering grid voltage, the SEs are rejected and only BSEs are detected.

Below a landing energy of 1.5 kV, the filtering grid **2** has the additional function of selecting the desired energy of the BSEs. The operator can select the threshold energy of inelastically scattered BSEs to enhance contrast and resolution.

The combination of InLens SE detector **3** and EsB detector **1** allows simultaneous imaging and mixing of a high contrast topography (SE) and a compositional contrast (BSE).

3.3.1.3 BSD Detector

Purpose The BSD detector is a pneumatically retractable backscattered electron detector. It is used for high efficiency material characterization even at low-kV applications. It has four separate diode quadrants S1 to S4 and an additional shadow segment S5 for strong topography contrast.




Fig. 20: BSD detector with four diode quadrants and a fifth shadow segment

NOTICE

Risk of property damage: Motorized specimen stage

Risk of damaging the detector when operating the motorized specimen stage.

- ▶ Retract the detector head completely after you have finished the work with the detector.

Info

Risk of malfunction: The diode segments are sensitive to the light that is used for illumination in TV mode (infrared and white).

When you use a diode detector, always make sure that the TV illumination is switched off. If the CCD Mode is set to Auto Detect, then the TV illumination is automatically switched off when a diode detector is used.

The BSD detector has applications mainly in materials analysis and in the life sciences.

Material analysis:

- Metallurgical sections
- Geological sections
- Complex materials
- Printed circuit boards
- Semiconductors
- Bond pads

Life sciences:

- Mineral deposits in plant structures
- Bone structures

The BSD detector is available either with one video output channel (BSD1) or with four video output channels (BSD4). With BSD4 you can collect four channels in parallel, which can be displayed in Quad mode. Each channel can have an arbitrary combination of segments. This way, you can collect compositional images, topographical images, and custom combinations all at once.

Four channels are especially useful in combination with the program 3DSM, which provides a topographical reconstruction of the specimen. The program 3DSM acquires four images of the specimen and acquires each of these images with a different quadrant of the detector. This results in four images from four different directions. With the BSD4 detector you can acquire these four images all at once. An algorithm then calculates the topographic surface.

The BSD detector has a relatively large central hole and therefore has the advantage that it does not limit the field of view of the SEM and does not influence the electron optical properties of the objective lens. A disadvantage is that especially at low kV a lot of backscattered electrons are lost in the central hole and cannot be detected. A solution for that is the aBSD detector.

Operating Principle On the surface of the specimen, some of the primary electrons are backscattered. The backscattered electrons then move towards the silicon segments of the BSD detector. If the energy of the backscattered electrons is high enough, then the electrons pass through the very thin dead layer of the diode and create electron-hole pairs in the silicon segments.

In each individual segment, the charge separation due to the electron-hole pairs is measured as a current, which is used as a signal for image generation. Only electrons that have a high enough energy can create electron-hole pairs and can contribute to image generation. Electrons that have a lower energy (e.g. secondary electrons) are not detected by the BSD detector.

The emission of backscattered electrons from a specimen is related to the atomic number of the involved material: Elements with high atomic numbers generate more backscattered electrons (i.e. the backscatter coefficient is higher). When imaging, regions that contain elements with higher atomic numbers appear brighter. Regions that contain elements with lower atomic numbers appear darker.

Since the detector has a limited speed, it is recommended to use scan speed 6 or higher (slower), especially at small magnifications. The lower the gain is, the faster is the detector.

3.3.1.4 aBSD Detector

Purpose The aBSD detector is a pneumatically retractable annular backscattered electron detector. It is used for high efficiency and angle selective material characterization even at low-kV applications. It has six separate diode segments, two inner concentric rings and four outer quadrants. The inner segments provide mostly material contrast whereas the four outer quadrants provide more topographical contrast.




Fig. 21: The aBSD detector with two inner concentric segments and four outer quadrants

NOTICE

Risk of property damage: Motorized specimen stage

Risk of damaging the detector when operating the motorized specimen stage.

- ▶ Retract the detector head completely after you have finished the work with the detector.

Info

Risk of malfunction: The diode segments are sensitive to the light that is used for illumination in TV mode (infrared and white).

When you use a diode detector, always make sure that the TV illumination is switched off. If the CCD Mode is set to Auto Detect, then the TV illumination is automatically switched off when a diode detector is used.

The BSD detector has applications mainly in materials analysis and in the life sciences.

Material analysis:

- Metallurgical sections
- Geological sections
- Complex materials
- Printed circuit boards
- Semiconductors
- Bond pads

Life sciences:

- Mineral deposits in plant structures
- Bone structures

The BSD detector is available either with one video output channel (BSD1) or with four video output channels (BSD4). With BSD4 you can collect four channels in parallel, which can be displayed in Quad mode. Each channel can have an arbitrary combination of segments. This way, you can collect compositional images, topographical images, and custom combinations all at once.

Four channels are especially useful in combination with the program 3DSM, which provides a topographical reconstruction of the specimen. The program 3DSM acquires four images of the specimen and acquires each of these images with a different quadrant of the detector. This results in four images from four different directions. With the BSD4 detector you can acquire these four images all at once. An algorithm then calculates the topographic surface.

The BSD detector has a relatively large central hole and therefore has the advantage that it does not limit the field of view of the SEM and does not influence the electron optical properties of the objective lens. A disadvantage is that especially at low kV a lot of backscattered electrons are lost in the central hole and cannot be detected. A solution for that is the aBSD detector.

Operating Principle

On the surface of the specimen, some of the primary electrons are backscattered. The backscattered electrons then move towards the silicon segments of the BSD detector. If the energy of the backscattered electrons is high enough, then the electrons pass through the very thin dead layer of the diode and create electron-hole pairs in the silicon segments.

In each individual segment, the charge separation due to the electron-hole pairs is measured as a current, which is used as a signal for image generation. Only electrons that have a high enough energy can create electron-hole pairs and can contribute to image generation. Electrons that have a lower energy (e.g. secondary electrons) are not detected by the BSD detector.

The emission of backscattered electrons from a specimen is related to the atomic number of the involved material: Elements with high atomic numbers generate more backscattered electrons (i.e. the backscatter coefficient is higher). When imaging, regions that contain elements with higher atomic numbers appear brighter. Regions that contain elements with lower atomic numbers appear darker.

Since the detector has a limited speed, it is recommended to use scan speed 6 or higher (slower), especially at small magnifications. The lower the gain is, the faster is the detector.




Fig. 22: Quad mode showing different contrast on different segments of the aBSD detector: Segment S1: Pure material contrast (top left), Segment S2: Material contrast and topography (top right), Segments S3 to S6: Pure topography (bottom left). For comparison: Surface contrast seen with the InLens SE detector (bottom right)

At small magnifications, the field of view is limited by the hole inside the aBSD detector:




Fig. 23: Central hole of the aBSD detector

The central hole is laser cut. As a result of this process, there is some silicon that beads up in areas near the hole, but this silicon is conducting. It is therefore not charging up and not disturbing your image quality.

NOTICE

Risk of property damage: Short working distance

When inserted, the aBSD detector is positioned directly underneath the objective lens. The lower edge of the detector is then located at a working distance of 1.5 mm. If you move the specimen to a working distance below 1.5 mm, then you damage the aBSD detector.

- ▶ Do not move the specimen to working distances below 2 mm.
- ▶ Be careful when you tilt the specimen.

3.3.1.5 aSTEM Detector

Purpose The optional aSTEM (annular Scanning Transmission Electron Microscopy) detector consists of an electron detector underneath an ultrathin specimen.

The aSTEM unit is equipped with diodes that are switched on or off in order to allow dark field and bright field imaging.




Fig. 24: aSTEM detector

- 1** Incident electron beam (primary electrons)
- 2** Thin specimen
- 3** aSTEM detector

Info

Risk of malfunction: The diode segments are sensitive to the light that is used for illumination in TV mode (infrared and white).

When you use a diode detector, always make sure that the TV illumination is switched off. If the CCD Mode is set to Auto Detect, then the TV illumination is automatically switched off when a diode detector is used.

Operating Principle The aSTEM detector is a pneumatically retractable multi-mode detector with a 12-stub specimen holder for bright-field and dark-field detection.

The diodes of the aSTEM detector collect transmitted electrons below the ultrathin specimen.

A special arrangement of the diodes allows a parallel detection of bright field (BF), oriented dark field (ODF), annular dark field (ADF), and high angle annular dark field (HAADF). There are one bright field segment (S1), two dark field segments (S2 and S3), one annular dark field segment (S4), and one high angle annular dark field segment (S5).

Two different models are available:

- The aSTEM4 allows detection of four channels in parallel.
- The aSTEM1 only collects one signal at a time.




Fig. 25: The aSTEM detector head and the arrangement of its diodes

Operating Modes

The most commonly used STEM modes are predefined and can be selected from the drop-down menu in the STEM Control panel.

The diagram in the STEM Control panel displays details concerning the diode status and the diode arrangement.

Each diode segment can have the following status:

- On: The signal is added to the total signal
- Inverted (Inv.): The signal is inverted and then added to the total signal
- Off: The signal is not used

NOTICE

Risk of property damage: Motorized specimen stage

Risk of damaging the detector when operating the motorized specimen stage.

- ▶ Retract the detector head completely after you have finished the work with the detector.

3.3.1.6 CL Detector

Purpose The Cathodoluminescence (CL) detector is an inclined detector that allows efficient visible or ultra-violet light collection. The CL detector is ideal for use in geology, mineralogy, and materials science applications where it can help in internal structural examination of rocks, ceramics, and semiconductors.




Fig. 26: Schematics of the CL detector

- 1 Preamplifier
- 2 Photomultiplier
- 3 Objective lens
- 4 Light guide with collector electrode
- 5 Specimen

Operating Principle The prerequisite for using this detector is that the specimen emits light when interacting with the primary electron beam. Differences in crystal structure or the presence of impurities in a cathodoluminescent material result in variations in the energy gap between the filled valence bands and the empty conduction bands, and consequently a change in the CL emission.

The light (photons) emitted by the specimen is collected by the CL detector and converted into a signal for imaging.

The CL detector is fully integrated into the automatic brightness and contrast control of the microscope and can be used simultaneously with any of the detectors without degrading their performance.

The detector can be used during energy-dispersive X-ray spectrometer (EDS) measurements and wavelength-dispersive spectrometer (WDS) measurements at any valid magnification.




Fig. 27: CL image of sandstone

3.3.2 Control Panel

Purpose The control panel allows direct access to 14 of the most frequently used functions. It integrates a full sized keyboard, 11 turning knobs, and 8 push buttons.

Position The control panel is placed on the work desk.




Fig. 28: Control panel

- 1** Stigmator X | Stigmator Y
Shapes the beam roundness by changing the stigmation deflectors.
- 2** Aperture X | Aperture Y
Adjusts the mid column shift and tilt deflectors for aligning the beam along the column axis.
- 3** Scan Rotate
Rotates the scanning pattern 360° continuously.
This turning knob has a push button function to deactivate the scan rotate function and reset the scan rotation to 0°.
- 4** Shift X | Shift Y
Shifts the scanned region of the specimen in the X and Y directions.
- 5** Brightness | Contrast

 - Brightness
Adjusts the image acquisition chain offset for the currently selected detector. Each configured detector stores its own brightness.
 - Contrast
Adjusts the gain of the currently selected detector.
- 6** Magnification | Reduced

 - Magnification
Adjusts the magnification of the system.
 - Reduced
Changes the scan field to a reduced area. The size of the area is determined by the current sub scan area settings.

7 Wobble

Sweeps the acceleration voltage. If the aperture is slightly misaligned, a shift in X and/or Y direction can be observed.

8 Freeze

Stops the scan and grabs one complete frame at the current imaging conditions.

9 Exchange | Resume

- Exchange

Starts the pre-defined macro for specimen exchange with the airlock.

- Resume

Starts the pre-defined macro to finish specimen exchange with the airlock.

10 Camera

Switches to chamber view.

11 Focus | Scan Speed +/-

- Focus

Changes the focal point of the column by adjusting the magnitude of the objective lens.

- Scan Speed +/-

Increases (+) or decreases (-) the scan speed by doubling or halving the beam dwell time with each click step.

4 Software Description

4.1 SmartSEM

4.1.1 SmartSEM User Interface

The SmartSEM software graphical user interface (GUI) allows you to monitor and operate most of the active components of the microscope.

The following screenshot indicates the main elements of the SmartSEM user interface:












Fig. 29: Screen layout of the user interface





- | | |
|--|---|
| <p>1 Title Bar</p> <p>Displays the name of the user interface and the logged-on user.</p> | <p>2 Menu Bar</p> <p>Enables you to access to SmartSEM features via sub-menus.</p> |
| <p>3 AVI Toolbar</p> <p>Contains the controls to set up, record, and play-back video sequences of scanned images.</p> | <p>4 Toolbar</p> <p>Provides quick access to SmartSEM tools.</p> |
| <p>5 Image Area with Data Zone</p> <p>Displays image information and acquisition parameters from the microscope.</p> | <p>6 FIB Toolbar</p> <p>Contains the controls to configure the FIB column.</p> |

- 7** **Thumbnails Panel**
Displays thumbnail views of the contents of the eight image buffers.
- 8** **Status Bar**
Displays the current machine state and contains the SEM Control Buttons.
- 9** **Data Zone**
Displays image information and acquisition parameters from the microscope.
- 10** **Annotation Bar**
Enables you to add information to the SEM image and provides several measurement functions.
- 11** **Panel Configuration Bar**
Enables you to choose the panels to be placed in the Docking Panel.
- 12** **Docking Panel**
Enables you to arrange frequently used SmartSEM panels for convenient access.
- 13** **Mini Bar**
Provides quick access to recently used dialogs and to the recipe management.

4.1.2 Graphical Control Elements

The following graphical control elements are used in the SmartSEM GUI.

| Screenshot | Control Element | Function |
|---|-----------------|--|
|  | Tab | Provides a group of graphical control elements. |
|  | Section | Forms a group of control elements with related functions. |
|  | Button | Enables you to start an action. |
|  | Checkbox | Enables you to activate or deactivate a function. |
|  | Drop-down list | Enables you to select the desired element. |
|  | Radio button | Enables you to activate the desired option. |
|  | Scroll bar | Enables you to adjust a value by moving the scroll bar or pressing the arrow button until the desired value is set. |
|  | Readout | Displays the status of a system entity. Enables you to select an action or a value by opening a dialog with an input field. |

| Screenshot | Control Element | Function |
|---|-----------------|--|
|  | Input field | Enables you to enter the desired value. |
|  | Progress bar | Displays the progress of an action. |
|  | Slider | Enables you to adjust the corresponding function. |
|  | Navigation box | Provides visual indication of the range and current value of one- and two-dimensional parameters such as Beam Shift or Stigmation. |

4.1.3 User Access Levels and User Privileges

The user access level defines which parameters are displayed for selection purposes, e.g. in the status window or annotation parameter selection.

SmartSEM distinguishes different user access levels. Depending on the user access level, different parameters are accessible. User profiles are defined by the administrator.

Access: Menu Bar > Tools > Administrator

| User Access Level | Description |
|-------------------|--|
| Novice | Only the items assigned to the novice category are accessible. These include most frequently used parameters. |
| Expert | Items assigned to the novice and expert category are accessible. These include parameters useful for advanced operators. |
| Service | All items are accessible, also including infrequently used items and calibrations. |

Additional to the user access levels there are user privileges which are part of the user profile:

| Checkbox | Privilege |
|------------------------|--|
| Calibration | Enables the user to perform instrument calibration operations. |
| Change Image Directory | Enables the user to change the location where all images are saved. |
| Change Toolbar | Enables the user to change the toolbar. |
| Change User Directory | Enables the user to change the location where all user specific parameters and configurations are saved. |
| Extractor | Enables the user to change the extractor voltage. |

| Checkbox | Privilege |
|---------------------|--|
| FIB Probe Alignment | Enables the user to adjust the probe currents. |
| Gun Align | Enables the user to modify the alignment of the electron beam. |
| Gun Off | Enables the user to switch off the field emission filament. |
| Mill Defaults | Enables the user to modify the default settings for FIB milling. |
| Stage Initialise | Enables the user to initialize the motorized stage. |
| Supervisor | <p>Enables the user to perform the following actions:</p> <ul style="list-style-type: none"> ▪ Start the Administrator, create and edit users ▪ Set User Max EHT ▪ Modify the filament current ▪ Set up, edit, and delete global stage coordinates ▪ Save common macros and toolbars ▪ Save common recipes ▪ Activate Partial Vent on Standby, Z Move on vent, Protect Z, Go to HV@Shutdown, EHT Off & Log Off and Leave Gun ON at Shutdown. ▪ Use the bakeout function ▪ Start the FIB filament heating. |

4.1.4 SmartSEM Program Suite

The SmartSEM Program Suite comprises the EM server, which implements the internal communication between control software and microscope hardware, plus several programs and utilities. The main purpose of the SmartSEM Program Suite is to access all necessary microscopy parameters and software features to capture SEM data and optimize image acquisition.

Access: Windows start menu > SmartSEM

| Program | Description |
|--------------------------|---|
| ChamberScope | Enables you to display the chamberscope image and the detector image at the same time. Option, requires particular hardware. |
| FTP Image Archiving | Enables you to transfer data via FTP. License: REMARCH |
| ReadMe | Contains important information on the currently installed version. |
| Release Notes | Contains an overview of all SmartSEM versions including new developments and specific details. |
| RemCon32 | Serial interface for remote operation via RS232, e.g. for EDX License: REMCON |
| Sample Holder Gallery | Enables you to inspect the dimensions of all possible specimen holders as well as to set the dimensions of the custom specimen holders. Enables you to activate the available specimen holders for SmartSEM. |
| SEM Drift Correction | Enables you to compensate for the drift of the specimen by using a reference image and by controlling the beam shift. License: DRIFT-CORR |
| Slideshow speed setting | Enables you to adjust the slideshow speed for the Windows Photo Viewer. |
| SmartSEM Administrator | Enables you to manage user profiles and configure instruments. |
| SmartSEM User Accounting | Enables you to record important information during individual working sessions, e.g. logon/logoff time, number of TIFF files exported etc. |
| SmartSEM User Interface | Main software application |

Access: Windows start menu > SmartSEM Service


| Program | Description |
|---------------------|--|
| Calibration Wizard | Service activities, for ZEISS service representatives only |
| Gun Monitor | Enables you to monitor important parameters of the microscope. |
| GUN Service | Service activities, for ZEISS service representatives only |
| Piezo Configurator | Service activities, for ZEISS service representatives only |
| Service Centre | Provides an overview of the state of the microscope. |
| Smart Stage Mapping | Service activities, for ZEISS service representatives only |

| Program | Description |
|-----------------------------|--|
| SmartBackup Tool | Enables you to back up configuration and calibration data. |
| Stage Administrator | Service activities, for ZEISS service representatives only |
| Upgrade Scangen Firmware | Service activities, for ZEISS service representatives only |
| Upgrade Server Database | Service activities, for ZEISS service representatives only |
| Merlin Alignment Wizard | Service activities, for ZEISS service representatives only |
| Merlin Database Wizard | Service activities, for ZEISS service representatives only |

4.2 SmartFIB

4.2.1 SmartFIB User Interface

The following screenshot indicates the main elements of the SmartFIB user interface:



1 Menu Bar

The menus on the Menu Bar contain basic commands that you need to work with SmartFIB. The functionality depends on the selected mode.

2 Standard Toolbar

Contains various buttons for quickly accessing a subset of the commands contained in the Menu Bar.

Additionally, there are various exclusive functions, such as previews and mode switching.

3 Working Area

In the Working Area, the following items are displayed:

- Live Mode: captured image and patterning elements
Editing of patterning elements is possible only in Live Mode.
- Sample Mode: Virtual Sample and writing positions
- EPD Mode: captured image and thickness map

4 Control Panel

Contains tabs for controlling the FIB process. The range of visible tabs depends on the selected mode.

5 Status Bar












- Left hand side: informs you about the current mouse position in the Working Area and indicates the selected tool, e.g. Vertex.
- Right hand side: enables you to control some exclusive properties of the Working Area.

5 Tools Toolbar

Contains tools for the Working Area. The range of visible tools depends on the selected mode.

4.2.2 Graphical Control Elements

The following graphical control elements are used in the SmartFIB GUI.

| Screenshot | Control Element | Function |
|---|--------------------|--|
|  | Tab | Provides a group of graphical control elements. |
|  | Section | Forms a group of control elements with related functions. |
|  | Button | Enables you to start an action. |
|  | Checkbox | Enables you to activate or deactivate a function. |
|  | Drop-down list | Enables you to select the desired element. |
|  | Radio button | Enables you to activate the desired option. |
|  | Scroll bar | Enables you to adjust a value by moving the scroll bar or pressing the arrow button until the desired value is set. |
|  | Readout | <ul style="list-style-type: none"> Displays the status of a system entity. Enables you to select an action or a value by opening a dialog with an input field. |
|  | Input field | Enables you to enter the desired value . |
|  | Slider | Enables you to adjust the corresponding function. |
|  | Expandable section | Enables you to expand or collapse the section as required. |

4.2.3 SmartFIB Program Suite

The SmartFIB program suite comprises two main programs: SmartFIB and Designer (requires the license CREATOR). Each of these programs has a different field of use.

The main purpose of the SmartFIB program suite is to transfer geometric elements to a specimen with the help of a particle beam. This process is referred to as "exposure".

Access: Windows Start Menu > Programs > SmartSEM

| Program | Description |
|----------|---|
| SmartFIB | <ul style="list-style-type: none"> ▪ Main tool for online/live work on the microscope ▪ Milling/etching/deposition of patterning elements ▪ Provides the following modes: <ul style="list-style-type: none"> – Live Mode: Mainly used for editing patterns and structures and creating recipes – Sample Mode: Mainly used for creating recurring/automated workflows – EPD Mode: Used for thickness determination of a TEM lamella |
| Designer | <ul style="list-style-type: none"> ▪ Offline creation of layouts: Arrangement of elements in the scanning area ▪ Interaction with Sample Mode: Used as a drawing tool for Sample Mode |

5 Installation

Unpacking, installation, and first start-up are carried out by an authorized ZEISS service representative.

An excerpt of the installation requirements can be found under *Installation Requirements* [[▶ 114](#)].

For more details refer to the document *Installation Requirements*.

6 Operation

6.1 Starting the System

6.1.1 Energizing the Microscope

WARNING

Risk of injury: Restart after emergency off

If the reason for the emergency off is not eliminated, it may be dangerous to restart the microscope.

- ▶ If the microscope has been de-energized due to an emergency, ensure that the reason for the emergency off does not exist anymore.
- ▶ Make sure it is safe to energize the microscope.

- Prerequisite**
- ✓ The microscope has been de-energized by turning the **MAIN** switch to the **OFF** position.
 - ✓ Alternatively, the microscope has been de-energized by pressing the **Emergency Off (EMO) button**.
1. Verify that the main shut-off valves for compressed air, nitrogen, and cooling water at the facility installation are operable.
Otherwise open and unlock the main shut-off valves.
 2. If the Emergency Off (EMO) option is installed, pull the **EMO** button to release it.
 3. Set the **MAIN** switch to the **ON** position.



4. If the Emergency Off (EMO) option is installed, press the green **S2** button to confirm the setting of the **MAIN** switch.



6.1.2 Starting the Microscope

Prerequisite ✓ The microscope is energized.

1. At the front of the plinth, press the **ON** button.
 - The **ON** button blinks green while the system is activated.
 - When all subsystems are fully activated, the **ON** button lights up green permanently.



6.1.3 Starting SmartSEM

1. Power up the computer and log in.
2. Start the SmartSEM user interface via the ZEISS SmartSEM icon on the desktop. Alternatively, select Windows start menu > SmartSEM > SmartSEM User Interface.
 - The EM Server opens while loading various drivers. The EM Server implements the internal communication between software and hardware of the microscope.
 - The EM Server Log On dialog is displayed.
3. Enter the user name and password.
4. Click **OK**.
 - The SmartSEM user interface opens and is ready to operate the microscope.

6.1.3.1 Calling up the Help

There are different ways to access topics in the Online Help.

| Function | Menu | Shortcut | Control Elements |
|-------------------|-----------------|---|--|
| Startup page | Help | F1 | – |
| Table of Contents | Help > Contents | Ctrl+F1 | – |
| Context-sensitive | – | <ul style="list-style-type: none"> ▪ Shift+F1 ▪ F1 on focus | Question mark icon in the main window and in modal dialogs |

Detailed information about using the help system is given in the Online Help directly.

6.1.3.2 Keyboard Shortcuts

The following keys are shortcut keys and have special meaning.

| Shortcut | Function |
|--------------|---|
| <F2> | Toggles Tool Bar on/off |
| <F2 + SHIFT> | Hysteresis removal |
| <F3> | Closes all windows except the Tool Bar and Status Bar |
| <F3 + SHIFT> | Toggles PC Plane ON/OFF |

| Shortcut | Function |
|--|--|
| <F4> | Step to next Magnification Table entry, or Undo Centre Feature Magnification |
| <F4 + CTRL> | Step to previous Magnification table entry |
| <F4 + SHIFT> | Exit from Magnification Table mode |
| <F5>, <F5 + SHIFT> <F6>, <F6 + SHIFT> <F7>, <F7 + SHIFT> <F8>, <F8 + SHIFT> | User defined macros |
| <F9> | Keys help (displays this information) |
| <F11>, <F11 + SHIFT> | User defined macros |
| <F12>, <F12 + SHIFT> | Aborts Stage Movement |
| <TAB> | Toggle coarse/fine |
| <CTRL + TAB> | Performs Centre Point |
| <CTRL + SHIFT + TAB> | Performs Centre Feature |
| <HOME> | Resets Beam Shift to zero |
| <SCROLL LOCK> | Toggles Freeze/Unfreeze |
| <PAUSE> | Causes currently executing macro to continue |
| <*> | Performs Find Image function |
| <CTRL + 2> | Loads Second Image Window from display |
| <CTRL + A> | Switches Annotation panel ON |
| <CTRL + B> | Display Toolbar View dialog |
| <CTRL + D> | Toggle Data Zone ON/OFF |
| <CTRL + E> | Calls the Export TIFF dialog |
| <CTRL + F> | Starts Auto Focus fine |
| <CTRL + SHIFT + F> | Starts Auto Focus coarse |
| <CTRL + G> | Switches Crossbeam SEM Control Panel ON |
| <CTRL + I> | Switches SEM Status Panel ON |
| <CTRL + M> | Switches to Annotation and inserts Point to Point Marker |
| <CTRL + ALT + M> | Enable/Disable the Movable Crosshairs Marker |
| <CTRL + ALT + F> | Enable/Disable Mouse Following for Movable Crosshairs Marker |

| Shortcut | Function |
|--------------------------|--|
| <CTRL + O> | Calls the Import TIFF dialog |
| <CTRL + P> | Performs the Print Image function |
| <CTRL + S> | Calls the Export TIFF dialog to save the image |
| <CTRL + ALT + S> | Performs Auto Astigmatism Correction |
| <CTRL + SHIFT + S> | Performs Auto Astigmatism Correction with Auto Focus |
| <CTRL + T> | Calls Text Annotation |
| <CTRL + V> | Displays the Vacuum status information |
| Keypad <+> | Faster Scan |
| Keypad <-> | Slower Scan |
| ARROW Keys | See Use of ARROW Keys |
| Image Buffer keys | See Image Buffer |
| <SHIFT> and double click | Performs Centre Point function |

6.2 Obtaining a First Image

Info

Mobile phones in the microscope room can cause image quality infringements and in worst case workflow interruptions.

Purpose This section describes basic procedures to obtain an image using the SE detector. To simplify the procedure, the description uses the Crossbeam SEM Control panel and status bar functions in the SmartSEM software.

Overview The procedure contains the following steps:

- *Preparing the Specimen Holder* [[▶ 69](#)]
- *Loading the Specimen Chamber* [[▶ 70](#)]
- *Locating the Specimen* [[▶ 72](#)]
- *Switching on the Gun* [[▶ 73](#)]
- *Switching on the EHT* [[▶ 74](#)]
- *Acquiring an Image* [[▶ 74](#)]
- *Optimizing the Image* [[▶ 76](#)]
- *Saving the Image* [[▶ 78](#)]

6.2.1 Preparing the Specimen Holder

| Parts and Tools | Designation | Part no. |
|-----------------|---|-------------------------------|
| | Allen key, 1.5 mm | Delivered with the microscope |
| | Stub | Delivered with the microscope |
| | Tweezers for specimen | Delivered with the microscope |
| | Specimen holder | Delivered with the microscope |
| | If necessary: carbon tape, conductive carbon, adhesive metal tape, or similar | – |
| | Appropriate specimen (with conducting properties, e.g. gold on carbon) | – |
| | Lint-free gloves | – |


NOTICE

Risk of property damage: Contamination caused by fingerprints


Contamination caused by fingerprints can lead to vacuum deterioration or prolonged pumping times.

- ▶ Always wear lint-free gloves when touching the specimen, specimen holder, or stage.


1. To attach a specimen to the stub, use conductive carbon, adhesive metal tape, carbon tape, or similar. Ensure that the specimen area that you want to analyze is in proper contact with the stub.



2. To insert the stub into the specimen holder, use tweezers.



3. To fix the stub to the specimen holder, tighten the location screw with the Allen key.



4. Note down which fix position is occupied by the specimen.

6.2.2 Loading the Specimen Chamber

Purpose

Info

If your microscope is equipped with the optional airlock, use the airlock for loading the specimen chamber. For more information refer to the respective instruction manual.

Overview The procedure contains the following steps:

- *Displaying the Crossbeam SEM Control Panel* [▶ 70]
- *Driving the Stage to a Low Position* [▶ 70]
- *Venting the Specimen Chamber* [▶ 70]
- *Mounting the Specimen Holder* [▶ 71]
- *Evacuating the Specimen Chamber* [▶ 72]

6.2.2.1 Displaying the Crossbeam SEM Control Panel

Prerequisite ✓ The SmartSEM user interface is started.

1. From the **Menu Bar**, select **Tools > Goto panel**.
 - ➔ The **Panel Configuration Bar** is displayed. It contains an alphabetical list of functions.
2. Double-right-click Crossbeam SEM Control.
 - ➔ The Crossbeam SEM Control panel is added to the docking panel.

6.2.2.2 Driving the Stage to a Low Position

Prerequisite ✓ The stage is initialized.

1. In the Crossbeam SEM Control panel, select the **Imaging** tab.
2. In the **Detector / Active Channel** section, select **TV** from the **Signal A** drop-down list.
 - ➔ The inside of the specimen chamber is visible in the **Image Area**.
3. In the Crossbeam SEM Control panel, select the **Stage** tab.
4. Activate the **Track Z** checkbox.
 - ➔ The current working distance (WD) is displayed in the **Data Zone**.
5. If the **Data Zone** is deactivated, activate it via **Menu Bar > View > Data Zone > Show Data Zone**.
6. Use the dual joystick to drive the specimen stage downwards to a low position.
 - NOTICE** **Risk of property damage: Observe the stage movement via camera to avoid crashing.**

6.2.2.3 Venting the Specimen Chamber

1. In the Crossbeam SEM Control panel, select the **Vacuum** tab.
2. Click **Vent**.
 - ➔ The **Vent** message box is displayed.
3. To start venting, click **Yes**.
 - ➔ The specimen chamber is purged with gaseous nitrogen.
 - ➔ If the Stage is not initialized, a system message is displayed. Refer to *Initializing the Stage* [▶ 99].

6.2.2.4 Mounting the Specimen Holder

Safety Information

⚠ WARNING**Suffocation hazard: Lack of oxygen**

Gaseous dry nitrogen is used to vent the specimen chamber during specimen exchange. Inhaling nitrogen may cause unconsciousness.

- ▶ During specimen exchange, keep the chamber door open as short as possible.
- ▶ Do not inhale the air from within the specimen chamber.
- ▶ Ensure that the area around the microscope is sufficiently vented.
- ▶ If you begin to experience symptoms of asphyxia (for example: rapid breathing, loss of mental alertness and/or muscular coordination, depression of sensations, emotional instability, fatigue) leave the room immediately and inform the facility's safety officer.

⚠ CAUTION**Risk of injury: Moving the specimen stage**

Fingers can be trapped by the moving specimen stage.

- ▶ Always close the chamber door before moving the specimen stage.
- ▶ To remove parts fallen into or near to the stage use a tool (e.g. tweezers) instead of your fingers.

⚠ CAUTION**Risk of injury: Closing the chamber door**

Fingers can be pinched when closing the chamber door.

- ▶ Use the door handle to close the chamber door.
- ▶ Ensure not to get your fingers caught in the chamber door gap.

NOTICE**Risk of property damage: Short working distance**

When opening the chamber door, the microscope or specimen can be damaged if the specimen stage is at a short working distance. If a BSD detector is inserted, it can be damaged as well.

- ▶ Always move the specimen stage to a long working distance before opening the chamber door.

NOTICE**Risk of property damage: Contamination caused by fingerprints**

Contamination caused by fingerprints can lead to vacuum deterioration or prolonged pumping times.

- ▶ Always wear lint-free gloves when touching the specimen, specimen holder, or stage.

1. Carefully open the chamber door.
2. If a specimen holder is mounted onto the specimen stage, remove it by sliding it out of the dovetail rails.
3. Mount the prepared specimen holder by sliding it into the dovetail rails.
Make sure that the dovetail is placed in the correct orientation so that the flat side of the dovetail of the specimen holder is flush with the milled edge of the specimen stage.

4. Carefully close the chamber door.
 - The specimen holder and the specimen inside the chamber are visible in the **Image Area**.

6.2.2.5 Evacuating the Specimen Chamber

1. In the Crossbeam SEM Control panel, select the **Vacuum** tab.
2. Click **Pump**.
 - Several vacuum status messages display the current vacuum levels.
 - As soon as the appropriate vacuum level is achieved, the vacuum status message **Vac Status = Ready** is displayed. This may take up to 5 minutes.

6.2.3 Locating the Specimen

Overview The procedure contains the following steps:

- *Positioning the Stub under the Electron Beam* [[▶ 72](#)]
- *Moving the Specimen to the Proper Height* [[▶ 73](#)]

6.2.3.1 Positioning the Stub under the Electron Beam

Safety Information

NOTICE

Risk of property damage: Driving the stage

While the stage is driven manually, there is a risk of damaging the objective lens and/or the specimen.

- ▶ Ensure not to hit the objective lens while driving the stage.
- ▶ Monitor the moving stage in TV mode.
- ▶ To stop the moving stage immediately, press **F12** or press the **Break** push button of the dual joystick panel.
- ▶ Manually lower the stage before you open the chamber door. Alternatively, activate the Z move on Vent checkbox in the Stage tab of the Crossbeam SEM Control panel.

1. In the Stage Navigation Bar, select Stage Sideview from the upper drop-down list and Stage Topview from the lower drop-down list.
INFO: To open the Stage Navigation Bar, navigate to View > Toolbars and activate Stage Navigation Bar (for Widescreen users). Alternatively, you can access the Stage Navigation Bar via Stage > Navigation.
2. Click Settings.
 - The Stage Navigation Settings dialog is displayed.
3. In the Stage Navigation Settings dialog, click Show Holder Gallery.
 - The Sample Holder Gallery dialog is displayed.
4. In the Sample Holder Gallery dialog, click the specimen holder you are using.
5. Activate the Is Available checkbox.
6. Close the Sample Holder Gallery dialog.
7. Close the Stage Navigation Settings dialog.
8. In the Stage Topview section of the Stage Navigation Bar, spot the stub with the specimen you want to observe.
9. To drive the stub directly under the electron beam, double-click the stub.

6.2.3.2 Moving the Specimen to the Proper Height

1. In the Stage Navigation Bar, drag the Zoom View slider to the right end, so that the schematics are zoomed in.
2. In the Crossbeam SEM Control panel, select the Imaging tab.
3. In the Detector / Active Channel section, select USB TV1 from the Signal A drop-down list.
 - The inside of the specimen chamber is visible in the Image Area.
4. Use the dual joystick to carefully move up the stage so that the stub you are using is in the center of the upper schematic.

NOTICE Risk of property damage: Observe the camera image in order not to crash into the pole piece.

- INFO: After loading the specimen, there is the possibility to use the Sample Type Selection Function to automatically set a number of key parameters. Refer to Using the Sample Type Selection Function.

6.2.4 Switching on the Gun

NOTICE

Risk of property damage: Schottky field emitter

If the Schottky field emitter is switched on and off too frequently or inappropriately, its lifetime is reduced.

- ▶ Avoid switching off the gun during the working week.
- ▶ Use Standby mode for the weekend or a break of up to a week.
- ▶ When using the Standby mode, activate the Partial Vent on Standby function.

Prerequisite ✓ The chamber and the gun head have been evacuated.

1. In the right part of the **Status Bar**, verify whether the gun is switched on or off.
 - If **Gun: ✓** or **All: ✓** is displayed, the gun is already switched on and you can skip the following steps.
 - If **Gun: ✗** is displayed, the gun is switched off. Follow the operating steps within this chapter.
2. In the Crossbeam SEM Control panel, select the **Vacuum** tab.
3. Verify that the **EHT Vac ready** readout is **EHT Vac ready = Yes**.
If not, the correct vacuum is not achieved. Check if the **Pump** procedure has been completed.
4. In the right part of the **Status Bar**, click **Gun: ✗**.
 - The pop-up menu for vacuum, gun, and EHT activation is displayed.
5. Click **Gun On**.
 - The gun runs up.
 - This may take up to 5 minutes.

6.2.5 Switching on the EHT

Purpose When you switch on the EHT, the gun starts emitting electrons.

- Prerequisite**
- ✓ The chamber and the gun head have been evacuated.
 - ✓ The gun has been switched on.
1. In the Crossbeam SEM Control panel, select the Control tab.
 2. Double-click the EHT Target readout.
 - The EHT Target window is displayed.
 3. In the input field, enter 10 and click OK.
 4. In the right part of the **Status Bar**, click **EHT: X**.
 - The pop-up menu for vacuum, gun, and EHT activation is displayed.
 5. Click EHT On.
 - The EHT runs up to 10 kV.
 - In the right part of the **Status Bar**, the vacuum, gun, and EHT status buttons merge to **All: ✓**.

6.2.6 Acquiring an Image

Purpose

Info

The following procedure describes the best way to quickly obtain an image without the control panel. You can also use the control panel to adjust aperture alignment, magnification/focus and brightness/contrast.

Overview The procedure contains the following steps:

- *Selecting the Column Mode and Probe Current [▶ 74]*
- *Selecting the SE Detector [▶ 74]*
- *Setting a Fast Scan Speed [▶ 75]*
- *Setting a Low Magnification [▶ 75]*
- *Setting a Long Working Distance [▶ 75]*
- *Adjusting Brightness and Contrast [▶ 75]*
- *Visualizing Details on the Specimen Surface [▶ 75]*

6.2.6.1 Selecting the Column Mode and Probe Current

1. In the Crossbeam SEM Control panel, select the **Control** tab.
2. In the **Column** section, click **Normal**.
 - The column mode is set to **Normal**.
3. In the **Beam** section, double-click the **I Probe** readout.
 - The **I Probe** window is displayed.
4. In the input field, enter 300.
 - The probe current is set to 300 pA.

6.2.6.2 Selecting the SE Detector

1. In the Crossbeam SEM Control panel, select the Imaging tab.

- In the **Detector / Active Channel** section, select **Signal A = SE2** from the **Signal A** drop-down list.
INFO: We recommend using the SE detector to obtain the first image. This detector provides a good signal-to-noise ratio even at long working distances.

6.2.6.3 Setting a Fast Scan Speed

- In the **Crossbeam SEM Control panel**, select the **Imaging tab**.
- From the **Scan Speed** drop-down list, select **Scan Speed = 1**.
INFO: The lower the scan speed number, the faster the electron beam scans across the specimen. **Scan Speed = 1** allows you to get an image quickly.

6.2.6.4 Setting a Low Magnification

- In the **Toolbar**, select the **Magnification+Focus/Auto Focus+Stig** icon.



→ The **Status Bar** displays the values for magnification and focus.

- In the **Status Bar**, click **Left: Mag =**.
→ The **Mag** window is displayed.
- In the **Mag** input field, enter **500**.
- Click **OK**.
→ The magnification is set to **Mag = 500 x**.


6.2.6.5 Setting a Long Working Distance

- In the **Status Bar**, click **Mid: WD =**.
→ The **WD** window is displayed.
- In the **WD** input field, enter **10**.
- Click **OK**.
→ The working distance is set to **WD = 10 mm**.

6.2.6.6 Adjusting Brightness and Contrast

- In the **Crossbeam SEM Control panel**, select the **Imaging tab**.
- In the **Detector / Active Channel** section, use the scroll bars to adjust brightness and contrast.

6.2.6.7 Visualizing Details on the Specimen Surface

- Select a detail on the specimen surface.
- Verify the **Magnification/Focus** function is activated.

- To adjust the magnification, hold down the left mouse button and drag the mouse within the **Image Area** in left/right direction.
→ The current magnification is indicated in the **Status Bar**.
- To adjust the focus, change the working distance. Hold down the mouse wheel and drag the mouse within the **Image Area** in left/right direction.
→ The current working distance is indicated in the **Status Bar**.
- Adjust contrast and brightness again.

6.2.7 Optimizing the Image

Purpose Once you have generated an initial image, you can adjust various parameters to optimize the image.




Info

The following procedure describes the best way to quickly optimize the image without the control panel. You can also use the control panel to adjust aperture alignment, magnification/focus and brightness/contrast.

Overview The procedure contains the following steps:

- *Adjusting the Magnification* [▶ 76]
- *Moving the Field of View at High Magnifications* [▶ 76]
- *Limiting the Scan Field* [▶ 77]
- *Aligning the Aperture* [▶ 77]
- *Selecting the Scan Speed* [▶ 77]
- *Correcting Astigmatism* [▶ 78]


6.2.7.1 Adjusting the Magnification

1. To switch to the Fine mode, in the Status Bar, click **Coarse** .
 - The **Coarse**  button changes to **Fine** .
2. Step by step, raise the magnification up to Mag 50,000 x and focus in between. To adjust the magnification and the focus, hold down the left mouse button or the mouse wheel, respectively, and drag the mouse within the Image Area.

6.2.7.2 Moving the Field of View at High Magnifications

Purpose If you want to move the field of view at high magnifications, use the Beam Shift function instead of moving the stage.

1. In the Crossbeam SEM Control panel, select the Control tab.
2. In the Alignment section, click Beam Offset.
3. To shift the beam, in the Beam Offset navigation box, use the scroll bars or the red marker.



6.2.7.3 Limiting the Scan Field

Prerequisite ✓ Adjusting the size and position of the small frame (reduced raster) requires the license REDUCED.

1. In the Toolbar, click the Reduced Raster/Apertures icon.



- A small scan frame is displayed. This frame defines the specimen area to be scanned by the electron beam.
 - The image outside the scan frame is frozen.
2. To change the position of the scan frame, click on the green border line and use the mouse to drag and drop the frame.
 3. To change the size of the scan frame, click on the small blue squares on the green border line and drag them to the desired size.
 4. Focus the image in the reduced raster.

6.2.7.4 Aligning the Aperture

1. In the Crossbeam SEM Control panel, select the Control tab.
2. In the Alignment section, click Focus Wobble.
INFO: Focus wobble is a function that sweeps the acceleration voltage. If the aperture is misaligned, a lateral and vertical shift can be observed.
→ The Focus Wobble window is displayed.
3. To adjust the wobble intensity, use the Wobble Amplitude scroll bar.
4. To accelerate the wobble speed, activate the Wobble Fast checkbox.
5. In the Control tab, click Aperture.
6. In the Aperture Align navigation box, use the scroll bars or the red marker to adjust the aperture alignment until there is no movement of the detail in X- and Y-direction.
INFO: The specimen detail should just be pulsating without shifting.
7. In the Focus Wobble window, click OFF to deactivate focus wobble.
→ The Focus Wobble window closes.
8. Refocus the image.

6.2.7.5 Selecting the Scan Speed

1. In the Toolbar, from the Faster/Slower drop-down list, select Scan Speed = 7.



Alternatively, in the Crossbeam SEM Control panel, select the Imaging tab, and from the Scan Speed drop-down list, select Scan Speed = 7.

- The scan speed is set to Scan Speed = 7.
2. Bring the image into focus.

6.2.7.6 Correcting Astigmatism

1. Ensure that the Reduced Raster function is active.
2. Select a detail (e.g. a mark or an edge) on the specimen surface.
Ensure that the selected detail is in the raster. You can move the stage or shift the beam for this purpose.
3. In the Crossbeam SEM Control panel, select the Control tab.
4. Click Stigmator.
5. In the Stigmation navigation box, use the scroll bars or the red marker to obtain the sharpest possible image.
INFO: The specimen detail should just be pulsating without shifting.
INFO: To obtain optimum results, alternately correct focus and astigmatism.
6. To deactivate the reduced raster, in the Toolbar, click the Reduced Raster/Apertures icon.


6.2.8 Saving the Image

1. In the Toolbar, click the Freeze:Unfreeze/Scanning icon.



→ A red dot at the right bottom of the Image Area indicates that the image is frozen.

2. From the Menu Bar, select File > Save Image.
→ The Export TIFF dialog is displayed.



3. To change the save path, click Change Directory.
→ A file explorer window is displayed.
4. To confirm the selected path, click Select Folder.
5. Enter the filename in the Filename input field.
6. Click Save <file name>.tif.
7. To continue imaging, click the Freeze:Unfreeze/Scanning icon.



6.3 Modifying Gun Parameters

6.3.1 Selecting the Gun Mode

Purpose The microscope can operate in different gun modes:

- Normal
- Imaging
- Analytic

Normal Gun Mode

In Normal gun mode, the temperature of the Schottky emitter (gun / filament) and the extraction voltage are optimized for a long lifetime of the Schottky emitter. Normal gun mode is suitable for most applications.

Imaging Gun Mode

In Imaging gun mode, the temperature of the Schottky emitter and the extraction voltage are reduced in comparison to the Normal gun mode. This leads to a reduction of the energy spread of the primary electrons. Overall, the probe current in Imaging gun mode is about half the probe current in Normal gun mode.

Imaging gun mode is especially useful at low kV or at high magnifications to reduce chromatic aberration and to achieve a better resolution. Switching from Normal gun mode to Imaging gun mode is useful for reducing the probe current without any need for beam adjustments.

Analytic Gun mode

In Analytic gun mode, the Schottky emitter runs at default temperature but the extraction voltage is increased in comparison to the Normal gun mode. This leads to a higher probe current. Overall, the probe current in Analytic gun mode is about twice the probe current in Normal gun mode.

The Analytic gun mode is especially useful for applications which require high intensities (e.g. WDX).

Info

After switching the gun mode, you can immediately work with the selected gun mode. For applications, which require a high probe current stability, wait 24 hours until a stability of 0.2 %/h is reached. It is recommended not to change the gun mode during quantitative specimen analysis.

1. From the Menu Bar, select Tools > User Preferences.
→ The User Preferences dialog is displayed.
2. Select User > Expert Gun Mode.
3. Click in the Value field and select Yes.
4. Close the User Preferences dialog.
5. In the Panel Configuration Bar, double-click Crossbeam SEM Control.
6. In the Crossbeam SEM Control panel, select the Control tab.
7. In the Beam section, click one of the following buttons:
To switch to Imaging gun mode, click Imaging.
To switch to Normal gun mode, click Normal.
To switch to Analytic gun mode, click Analytic.
INFO: For maximum probe current, also check the alignment of the gun.
INFO: If you use the Analytic gun mode, then the lifetime of the emitter is reduced.

6.3.2 Setting the Probe Current

Purpose With the GEMINI II column, you can set a lower probe current to analyze surface details at a high resolution or higher probe currents for analytical purposes, e.g. to analyze the material of the specimen.

Info

The maximum achievable probe current depends on the currently selected EHT and the installed aperture configuration.

1. In the Crossbeam SEM Control panel, select the Control tab.
2. Double-click the I Probe readout.
 - The I Probe window is displayed.
3. In the input field, enter the desired value.

6.3.3 Measuring the Probe Current

Purpose Measuring the probe current using the Faraday cup ensures that the current displayed in the software equals the incident probe current. The Faraday cup consists of a strongly absorbing material with a cavity covered by a small aperture. If the beam is focused in this cavity, no secondary electrons and no backscattered electrons leave the Faraday cup.

Parts and Tools

| Designation | Part no. |
|-------------|-----------------|
| Faraday cup | 348342-8055-000 |

1. Load the Faraday cup into the specimen chamber.
2. Evacuate the specimen chamber.
3. Switch on the gun.
4. Switch on the EHT.
5. From the Panel Configuration Bar, select Specimen Current Monitor.
 - The Specimen Current Monitor window is displayed.
6. Activate the Stage Bias checkbox.
 - This activates the touch alarm that helps to avoid collisions of the stage.
7. Move the stage to the position of the Faraday cup.
8. Acquire an image of the Faraday cup.
9. Activate the Spot checkbox.
 - Green crosshairs are displayed on the image. The crosshairs indicate the position of the beam spot.
10. Grab the crosshairs and move them into the hole of the Faraday cup.
 - The probe current is measured continuously.
 - The measured probe current is displayed in the Specimen I readout.

6.3.4 Changing the Extractor Voltage

Purpose The Extractor voltage is preset by the factory or by the ZEISS service representative. Within certain limits, the operator may carefully increase the extractor voltage in order to optimize the probe current for particular applications.

Info

Use a Faraday cup to measure the probe current when changing the extractor voltage.

Info

The newly set extractor value is only valid for the current work session. After a restart of the SmartSEM software, the microscope restores the nominal extractor voltage.

Info

Reducing the extractor voltage may impair the performance and resolution of the microscope.

- Avoid reducing the extractor voltage.
- If at all, reduce the extractor voltage only for a short time (1–2 h) and by maximum 500 V.

Prerequisite ✓ The user privilege Extractor is required to change the extractor voltage.

1. From the Menu Bar, select Beam > Gun Setup.
 - ➔ The Gun Service dialog is displayed.
2. To increase the extractor voltage, double-click the Extractor V Target readout.
 - ➔ The Extractor V Target window is displayed.
3. Enter a higher value.
4. Click OK.

6.4 Working with Different Aperture Configurations and Beam Modes

6.4.1 Determining the Installed Aperture Configuration | Gemini II Column

Purpose The achievable maximum probe current depends on the type of installed anode aperture. The type of aperture installed on the microscope can be determined via SmartSEM.

For the Gemini II column, two different column configurations are available:

- 40 nA high resolution configuration

| Anode aperture diameter | Probe current | Typical application |
|-------------------------|----------------|---------------------|
| 55 µm* | 10 pA to 40 nA | High resolution |

- 100 nA high current configuration

| Anode aperture diameter | Probe current | Typical application |
|-------------------------|-----------------|--|
| 90 µm* | 10 pA to 100 nA | Combined high resolution and analytical investigations |

* Calibration value: deviation of 10 % possible

Info

If you wish to change the installed configuration of your microscope, contact your local ZEISS service representative.

1. From the Menu Bar, select View > SEM Status.
 - The SmartSEM Status dialog is displayed.
2. In the Select tab, click Anode Aperture Diameter.
3. Select the Display tab.
 - The parameter Aperture Size is displayed.

6.4.2 Selecting the Column Mode

Purpose With the Gemini column, different column modes are available.

The Normal mode is the standard imaging and analytical mode with a high flexibility for different applications. Overview mode provides a large field of view.

1. In the Crossbeam SEM Control panel, select the Control tab.
2. In the Column section, select between:
 - Normal and Overview
3. If required, adjust the Depth of Field slider.


6.5 Finding Appropriate Detector Settings**6.5.1 Selecting a Detector**

You need to select an appropriate detector depending on the application and the pressure mode. In addition to the standard SE detector, several optional detectors are available.

For information on special set-up procedures for the detectors, see:

- *Setting up the InLens SE Detector* [▶ 83]
- *Setting up the SE Detector* [▶ 83]
- *Setting up the SESI Detector* [▶ 85]
- *Setting up the aBSD/BSD Detector* [▶ 86]
- *Setting up the EsB Detector* [▶ 84]
- *Setting up the aSTEM/STEM Detector* [▶ 88]
- *Setting up the CL Detector* [▶ 91]

1. Select the Imaging tab of the Crossbeam SEM Control panel.
2. Select the detector from the Signal A drop-down list.



6.5.2 Setting up the InLens SE Detector

Purpose The InLens SE detector collects the SE signal, acquiring mainly information about surface topography.




Fig. 30: Silver nanoparticles embedded in zeolite, imaged at 1.5 kV.

The following settings are recommended for the InLens SE detector:

| EHT | Typical WD | Recommended WD |
|---------------|------------|--|
| 20 V – 10 kV | 0–5 mm | Short working distances are preferable for good detection efficiency |
| 10 kV – 20 kV | 2–5 mm | |

Info

Avoid strong specimen tilting for the InLens SE detector.

1. In the Crossbeam SEM Control panel, select the Imaging tab.
2. From the Signal A drop-down list, select InLens.
3. Adjust the EHT and the working distance (WD) according to the suggestions in the table in order to optimize the image.

6.5.3 Setting up the SE Detector

Purpose The SE detector collects the SE signal, highlighting the topography of the specimen.




Fig. 31: An eledone tentacle

The following settings are recommended for the SE detector:

| EHT | Typical WD | Collector Voltage |
|--------------|------------|--|
| 500 V – 5 kV | 2–8 mm | <ul style="list-style-type: none"> Adjustable from –250 V to +400 V Standard applications: +300 V <p>At a high magnification, you can optimize the image by varying the collector voltage.</p> <ul style="list-style-type: none"> Pseudo-backscattered (BSE) image: –250 V to –50 V <p>This produces a topographic image of the sample with no material contrast.</p> |
| 5 kV – 30 kV | min. 6 mm | |

1. In the Crossbeam SEM Control panel, select the Imaging tab.
2. From the Signal A drop-down list, select SE2.
3. Adjust the EHT, working distance (WD), and collector voltage according to the suggestions in the table in order to optimize the image.

6.5.4 Setting up the EsB Detector

Purpose The EsB detector can be used to collect the backscattered electrons (BSE) signal. The BSE signal contains information about the material contrast. In the final image, heavy elements are represented by brighter pixels and light elements are represented by darker pixels.

By adjusting the filtering grid, energy-selected BSE images can be obtained. If the filtering grid voltage is set to 0, SE and BSE mixed images can be acquired.




Fig. 32: Silver nanoparticles embedded in zeolite, imaged at 1.5 kV

The following settings are recommended for the EsB detector:

| EHT | Typical WD | Filtering Grid |
|---------------|------------|---|
| 500 V – 10 kV | 0–5 mm | EsB Grid > 400 V to filter out the SE signal |
| 20 V – 500 V | 0–3 mm | EsB Grid = 0 V for use as an additional SE detector |

1. In the Crossbeam SEM Control panel, select the Imaging tab.
2. From the Signal A drop-down list, select EsB.
3. Adjust the EHT, working distance (WD), and EsB Grid according to the suggestions in the table in order to optimize the image.

6.5.5 Setting up the SESI Detector

Purpose The SESI detector is optionally available and replaces the chamber SE detector.

The SESI detector enables you to acquire both secondary electron images and FIB secondary ion images.




Fig. 33: Intergranular corrosion in an nickel based superalloy

The following settings are recommended for the SESI detector:

- Settings when working in SE mode, secondary electron imaging, FIB mode = SEM:

| EHT | Typical WD | Collector Voltage |
|---------------|---------------------------|--|
| 100 V – 30 kV | 2–12 mm Typically 5 mm | <ul style="list-style-type: none"> ▪ Adjustable from 0 V to +1500 V ▪ Best detection: +300 V to +400 V |

- Settings when working in SE mode, secondary electron imaging, FIB mode = FIB:

| EHT | Typical WD | Collector Voltage |
|--------------|-------------------|--|
| 2 kV – 30 kV | Coincidence point | <ul style="list-style-type: none"> ▪ Adjustable from 0 V to +1500 V ▪ Best detection: +300 V to +400 V |

- Settings when working in ion mode, secondary ions imaging, FIB mode = FIB:

| EHT | Typical WD | Collector Voltage |
|--------------|-------------------|--|
| 2 kV – 30 kV | Coincidence point | <ul style="list-style-type: none"> ▪ Adjustable from –4 kV to +0 kV ▪ Best detection: Around –4 kV |

1. In the FIB Toolbar, from the Imaging Mode drop-down list, select an imaging mode, e.g. FIB mode = SEM.
2. In the Crossbeam SEM Control panel, select the Imaging tab.
3. From the Signal A drop-down list, select SESI.
 - ➔ By default, secondary electrons are detected.
4. In order to detect secondary ions, in the Imaging tab, activate the SESI Mode checkbox.
5. Adjust the EHT, working distance (WD), and collector voltage according to the suggestions in the table in order to optimize the image.

6.5.6 Setting up the aBSD/BSD Detector

Purpose The aBSD detector is a pneumatically retractable backscattered electron detector which is inserted below the objective lens and is used for high efficiency and angle selective (aBSD) material characterization.





Fig. 34: Quad mode showing different contrast on different segments of the aBSD detector: Segment S1: Pure material contrast (top left), Segment S2: Material contrast and topography (top right), Segments S3 to S6: Pure topography (bottom left). For comparison: Surface contrast seen with the InLens SE detector (bottom right)

The following settings are recommended for the aBSD detector:

| EHT | Typical WD | Detector Settings |
|-----------|------------|--|
| 0.5–30 kV | 5–12 mm | Compositional mode for obtaining an image that is high in atomic number contrast Use topographic mode for showing surface details Use individual settings for channeling contrast Make use of different concentric rings of the aBSD detector to get angular resolved BSE images Use high or very high detector gain for low accelerating voltage and/or low beam current The EHT should be bigger than 2 kV to achieve a significant detection efficiency. |

- From the Panel Configuration Bar, select BSD Control. The BSD Control panel enables you to change the polarity of the segments, select BSD modes, and set the BSD gain.

→ The BSD Control panel is displayed.



2. In the BSD Control panel, click BSD in to insert the aBSD detector.
 - The stage is lowered by 20 mm to give space for the detector to be inserted.
 - The detector is inserted.
3. In the Crossbeam SEM Control panel, select the Imaging tab.
4. From the Signal A drop-down list, select aBSD1 or aBSD4 A depending on your detector type and configuration.

INFO: A refers to the video channel. B, C, and D are also available if you have the four-channel version.
5. Adjust the EHT and the working distance (WD) according to the suggestions in the table in order to optimize the image.
6. In the BSD Control panel, click a quadrant symbol to toggle its status between on (white), inverted (black), and off (gray).
7. To choose the respective mode, click BSD: COMPO or BSD: TOPO.

The default setting for BSD is BSD: COMPO. All six segments (S1–S6, aBSD) or four segments (S1–S4, BSD) are set to **on** and an image that is high in atomic number contrast is obtained.

The default mode for BSD: TOPO is S3 on, S4 on, S5 inv and S6 inv (aBSD) and S1 on, S2 on, S3 inv and S4 inv (BSD).
8. If you want to change the default setting to BSD: TOPO, click BSD: Set TOPO.
9. From the BSD Gain drop-down list, select Low, Medium, High, or Very High.

INFO: Since the detector has a limited speed, it is recommended to use scan speed 6 or higher (slower), especially at small magnifications. The lower the gain is, the faster is the detector.

6.5.7 Setting up the aSTEM/STEM Detector

Purpose The aSTEM/STEM detector is used for compositional imaging or topographical imaging of ultra-thin specimens. The aSTEM/STEM detector is available either with one video output channel (aSTEM1) or with four video output channels (aSTEM4). The aSTEM detector is optionally available.

The aSTEM/STEM detector is equipped with several separate diode segments. The signals of the segments can individually be added to or subtracted from the output signal in order to allow different STEM imaging modes, e.g. bright field (BF) or oriented dark field (ODF). The most commonly used STEM imaging modes are predefined and can be selected from a drop-down menu in the STEM Control panel.






Fig. 35: Asbestos, oriented darkfield (ODF) image captured with an aSTEM4 detector (left), ODF image showing real information about bending and lattice defects within each fiber (right).

| STEM Control | STEM imaging mode | Output signal |
|--------------|--------------------------------|---|
| | BF | S1 |
| | DF | S2 + S3 |
| | ODF | S2 - S3 |
| | ODF90° | S3 - S2 |
| | ODF + BF | S1 + S2 - S3 |
| | ODF90° + BF | S1 - S2 + S3 |
| | BF + DF | All on |
| | User defined name | Custom combination of segments |
| | Active segments with +/- signs | Custom combination of segments. "+" on, "-" inverted |

| STEM Control | STEM imaging mode | Output signal |
|--|--------------------------------|---|
|  | BF | S1 |
| | DF | S2 + S3 |
| | ODF | S2 - S3 |
| | ODF90° | S3 - S2 |
| | ODF + BF | S1 + S2 - S3 |
| | ODF90° + BF | S1 - S2 + S3 |
| | BF + DF | S1 + S2 + S3 |
| | ADF | S2 + S3 + S4 |
| | HAADF Middle | S4 |
| | HAADF outer | S5 |
| | HAADF | S4 + S5 |
| | User defined name | Custom combination of segments |
| | Active segments with +/- signs | Custom combination of segments. "+" on, "-" inverted |

| STEM Control | STEM imaging mode | Output signal |
|--------------|--------------------------------|--|
| | Off | - |
| | BF | S1 |
| | ADF | S2 + S3 + S4 + S5 |
| | ODF 0° | S2 - S3 + S4 - S5 |
| | ODF 90° | S3 - S2 - S4 + S5 |
| | HAADF | S6 |
| | ADF + HAADF | S2 + S3 + S4 + S5 + S6 |
| | DF 0° | S2 |
| | DF 90° | S3 |
| | DF 180° | S4 |
| | DF 270° | S5 |
| | User defined name | Custom combination of segments |
| | Active segments with +/- signs | Custom combination of segments. "+" on, "-" inverted |

- In the Crossbeam SEM Control panel, select the Imaging tab.
- For displaying several channels of a STEM/STEM detector simultaneously, from the Menu Bar, select Scanning > Quad Mode.
 - The Image Area is divided into 4 zones. To select a detector for a zone, click in the zone. An anchor symbol is displayed in the selected zone.
- From the Signal A drop-down list, select a STEM detector, e.g. aSTEM1.
- In the Panel Configuration Bar, double-click STEM Control.
 - To open the Panel Configuration Bar, from the Menu Bar, select Tools > Goto Panel.
 - The STEM Control dialog is displayed. In the upper section, the STEM Control dialog displays the status of the diode segments. The status is either on (white), inverted (black), or off (gray).
- Either select a STEM imaging mode from the STEM Seg. Mode/STEM 6 Seg. Mode drop-down list or click a custom selection of diode segments to toggle its status between on, inverted, and off.
- Set the Gain range. Select between Low, Medium, High, or Very High.
- For displaying several channels of a STEM/STEM detector simultaneously, repeat steps 3 to 6 for the other display zones.

6.5.8 Setting up the CL Detector

Purpose The CL detector is optionally available.

The CL detector collects visible or ultraviolet light and is especially useful for internal structural examinations of rocks, ceramics, and semiconductors.




Fig. 36: Zircon grains

The following settings are recommended for the CL detector:

| EHT | Typical WD |
|--------------|------------------------|
| 5 kV – 30 kV | 6–10 mm (min. 4 mm) |

1. In the Crossbeam SEM Control panel, select the Imaging tab.
2. From the Signal A drop-down list, select CL.
3. Adjust the EHT and the working distance (WD) according to the suggestions in the table in order to optimize the image.

6.6 Shutting down the System

6.6.1 Finishing the Work Session

NOTICE

Risk of property damage: Schottky field emitter

If the Schottky field emitter is switched on and off too frequently or inappropriately, its lifetime is reduced.

- ▶ Avoid switching off the gun during the working week.
- ▶ Use Standby mode for the weekend or a break of up to a week.
- ▶ When using the Standby mode, activate the Partial Vent on Standby function.

Prerequisite ✓ The microscope is in ON mode.

1. In the Crossbeam SEM Control panel, select the Gun tab.
2. Activate the EHT Off @ Log Off checkbox.
3. Activate the Leave Gun On at Shutdown checkbox.
4. Only when interrupting work for longer periods between 2 and 7 days: In the Vacuum tab activate the Partial Vent on Standby checkbox.
 - ➔ This maintains the gun vacuum, and switches off and protects the turbo pump.

5. Exit the SmartSEM user interface. Refer to *Exiting the SmartSEM User Interface* [▶ 92].
6. Close all programs and software.
7. Select Windows start button > Power icon > Shut down.
8. At the front of the plinth, press the **Standby** button.
 - The microscope switches to Standby mode.

6.6.1.1 Exiting the SmartSEM User Interface

1. From the Menu Bar, select File > Exit.
 - A system message is displayed.
2. Click Yes.
3. Close the EM Server.
 - A system message is displayed.
4. Click Yes.

6.6.2 De-energizing the Microscope

Prerequisite ✓ The microscope is in Standby mode, see *Finishing the Work Session* [▶ 91]

1. To switch to Off mode, press the **Off** button at the front of the plinth.
 - The yellow **Standby** button blinks during shut down.
 - When all subsystems are fully deactivated, the **Off** button lights up red permanently.
 - Computer, electronic components, and vacuum system are switched off.
 - The electron optical column is partially vented.
 - A 24 V auxiliary voltage is still present to restart the microscope.
2. To cut off the compressed air supply, close the compressed air main shut-off valve.
3. To cut off the nitrogen supply, close the nitrogen main shut-off valve.
4. To cut off the cooling water supply, close the main shut-off valve for water.
5. Set the **MAIN** switch to the Off position.
 - All power is cut off from the microscope.
6. To disconnect the power connection at the facility installation, remove the power plug (male connector) from the socket (female connector).
7. To secure the microscope against accidental re-activation, refer to *Performing a Lockout/Tagout* [▶ 93].

See also

- 📖 [Finishing the Work Session](#) [▶ 91]
- 📖 [Performing a Lockout/Tagout](#) [▶ 93]

6.6.3 Performing a Lockout/Tagout

The Lockout/Tagout (LOTO) procedure defines the process by which the microscope is isolated from potential hazards associated with the unexpected release of hazardous energy. This procedure needs to be performed when service, maintenance work, or unscheduled repairs can cause exposure to any form of hazardous energy.

Info

The customer is responsible for providing instructions on how to operate the energy-isolating devices properly.

Prerequisite ✓ The microscope has been de-energized.

1. Verify that the microscope has been de-energized properly.
2. To secure the microscope against accidental re-activation, lock the energy-isolating devices.

6.7 Performing an Emergency Shutdown

NOTICE

Risk of property damage: Components in the high voltage circuitry

When the microscope, especially the gun, is fully on, an abrupt shutdown of all electrical supplies may damage some components in the high voltage circuitry, mainly the cathode.

- ▶ Use the emergency off only in an emergency situation with personnel injury.

1. Depending on whether the EMO option is installed on your microscope or not, operate the respective control element.
If the EMO option is installed, press the **EMO** button on the top of the plinth.
If the EMO option is not installed, set the **MAIN** switch at the rear side of the plinth to the OFF position.

7 Maintenance and Repair

7.1 Maintenance Work

The preventive maintenance is performed by the ZEISS service representative and includes the following items:

- Inspection
- Preventive maintenance work
- Change of consumables and chemicals
- Equipment test
- Verification run

Info

The maintenance work is accomplished according to standardized maintenance plans and is recorded by the ZEISS service representative.

7.2 Maintenance Intervals

The maintenance intervals depend on the period of application of the device:

- 24 hours, 7 days a week: semiannually
- 8 hours, 5 days a week: annually

Info

Keep track of maintenance work and contact the ZEISS service representative in time.

A list of ZEISS locations and authorized service partners can be found at:

<http://www.zeiss.com/microscopy>

7.3 Change of Consumables and Chemicals

The change of consumables and chemicals has to be performed by a ZEISS service representative at mandatory intervals.

The times scheduled are designed for the maximum equipment performance level (i.e. 24 h per day of permanent operation).

| Interval | Component/Part |
|---|---|
| Every 6000 h ¹ (filament on) | Field emission gun (filament) ¹ |
| Yearly or as required (and after exchange of field emission gun) | <ul style="list-style-type: none"> ▪ Multihole/singlehole aperture ▪ Anode aperture ▪ Extractor aperture ▪ Anode aluminum seal ▪ Copper seal at gun head |
| Yearly or as required | Tip seal of the pre-vacuum pump |
| Yearly performance check | <ul style="list-style-type: none"> ▪ SE detector ▪ InLens SE detector |

| Interval | Component/Part |
|-------------------------------|---|
| | <ul style="list-style-type: none"> ▪ EsB detector |
| As required, approx. 3000 µAh | <ul style="list-style-type: none"> ▪ Ion source ▪ FIB apertures |
| If used up | Precursors |

¹ *There is no warranty on field emission guns; manufacturers do not guarantee any lifetime.*

Tab. 3: Schedule for the change of consumables

8 Troubleshooting

8.1 Overview

The following table provides hints for solving common problems. If you cannot solve the problem or if you are unsure about a certain technical difficulty, contact your local ZEISS service representative.

| Keyword | Symptom | Cause | Remedy |
|---------|---|--|--|
| Drift | Specimen seems to be moving. | <ul style="list-style-type: none"> Charging effects. Nonconductive specimen. | <ul style="list-style-type: none"> Ensure proper conduction of the specimen. Optimize specimen preparation. Apply a charge compensation method. |
| | | Stub not correctly fixed by screw. | Fix the stub correctly. |
| EHT | EHT cannot be switched on. | CAN communication has failed. | Refer to <i>Checking the CAN Communication</i> [▶ 99]. |
| | The workstation has crashed. | CAN communication has failed. | Refer to <i>Checking the CAN Communication</i> [▶ 99]. |
| FIB | Emission is unstable. | Ion source needs to be regenerated. | Refer to <i>Regenerating the Ion Source by Heating</i> [▶ 106]. |
| | | Ion source may be used up (> 3000 µAh). | Contact the ZEISS service representative to have the ion source replaced. |
| FIB | No emission. | Ion source needs to be regenerated. | Refer to <i>Regenerating the Ion Source by Heating</i> [▶ 106]. |
| | | Ion source may be used up. | Contact the ZEISS service representative to have the ion source replaced. |
| FIB | No image. | Aperture position is not correct. | Initialize the FIB aperture in the FIB tab of the FIB Control panel. Align the aperture in the Align tab of the FIB Control panel. |
| | | Column conditions have deteriorated, charging effects occur. | Start a purging process in the Options tab of the FIB Control panel. |
| | | Gun valve is closed. | Open the gun valve in the FIB tab of the FIB Control panel. |
| FIB | Emission current is too high (e.g. 4 or 5 µA) for more than 30 min. | Ion source needs to be regenerated. | Refer to <i>Regenerating the Ion Source by Heating</i> [▶ 106]. |

| Keyword | Symptom | Cause | Remedy |
|---------------|---|---|--|
| | | Automatic regulation of the emission current is deactivated. | Activate the Regulate checkbox in the FIB tab of the FIB Control panel. |
| | | Extractor target is not set to an appropriate value. | Change the extractor target. Refer to <i>Manually Finding a New Extractor Value</i> [▶ 107]. Open the Gun Monitor and record the parameters FIB Suppressor Target and FIB Emission I and send it to your ZEISS service representative. |
| FIB | Suppressor voltage is approaching 0 V or 2 kV. | Ion source needs to be re-generated. | Refer to <i>Regenerating the Ion Source by Heating</i> [▶ 106]. |
| FIB | Misalignment of structures when milling with different FIB probe currents. | Beam shift for the probe currents is not adjusted. | Adjust the beam shift correction for the relevant FIB probe currents. |
| Image quality | Image quality gets worse, but there is no change in total emission current. | Field emission gun has been damaged due to arcing. | Contact your local ZEISS service representative to have the field emission gun replaced. |
| | Image is noisy and noise reduction methods do not help. | Field emission gun is used up. | Contact your local ZEISS service representative to have the field emission gun replaced. |
| | Image is bad at low EHT (e.g. 1 kV) | Working distance is too long. | Reduce the working distance to a maximum of 7 mm. |
| InLens image | InLens image is noisy. | Working distance is too long. | Reduce the working distance. |
| | No InLens image can be obtained. | EHT exceeds 20 kV. | Reduce EHT to a maximum of 20 kV. |
| Microscope | Microscope is dead. | Circuit breaker is tripped (lower position). | Refer to <i>Checking the Position of the Circuit Breakers</i> [▶ 104]. |
| PC | Stored position of the specimen stage cannot be approached correctly. | PC has crashed. | Restart the PC. |
| | | Stage needs to be driven to a well-defined position. | Refer to <i>Initializing the Stage</i> [▶ 99]. |
| SEM Gun | Gun is switched off automatically. | Gun has been switched off automatically for safety reasons since gun vacuum is worse than 2×10^{-8} mbar | Refer to <i>Baking out the Gun Head</i> [▶ 102]. |

| Keyword | Symptom | Cause | Remedy |
|--|---|---|---|
| SE image | SE image is noisy. | Scintillator is used up. | Contact your local ZEISS service representative to have the scintillator replaced. |
| Specimen current | Specimen current is low. | Field emission gun is used up. | Contact your local ZEISS service representative to have the field emission gun replaced. |
| | | Working distance is too short. | Enlarge working distance to about 5 mm or more. |
| Specimen stage | Stage does not move. | Stage needs to be initialized. | Refer to <i>Initializing the Stage</i> [▶ 99]. |
| | Stage does not move accurately. | Stage needs to be initialized. | Refer to <i>Initializing the Stage</i> [▶ 99]. |
| | Stored position of the specimen stage cannot be approached correctly. | Absolute stage movement is required. Stage needs to be driven to a well-defined position. | Refer to <i>Initializing the Stage</i> [▶ 99]. |
| Stage/Joystick | Under TV control, the direction of dual joystick movement and direction of stage movement seem to be different. | TV joystick angle does not fit for the selected CCD camera. | Refer to <i>Changing the Joystick TV Angle</i> [▶ 100]. |
| | Stage cannot be moved by using the joystick. | Joystick Disable checkbox is activated. | Deactivate the checkbox in the Stage tab of the Crossbeam SEM Control panel. |
| Temperature, water flow | Error message Stage Board too hot (or similar) is displayed. | Flow of cooling water is not OK. | Refer to <i>Checking the Water Flow and Temperature</i> [▶ 102]. |
| Touch alarm | Touch alarm message is displayed. | Specimen or specimen holder has touched objective or wall of the specimen chamber. | Refer to <i>Resetting the Touch Alarm</i> [▶ 101]. |
| Vacuum | Vac ready = OK is not displayed after specimen exchange. | System vacuum is bad due to a vacuum leak at the chamber door. | Check the chamber door seal for cleanliness. If required, refer to <i>Replacing the Chamber Door Seal</i> [▶ 101]. |
| | Vac ready = OK is displayed very late after specimen exchange. | Gas ballast at rotary pump or scroll pump is activated. | Deactivate gas ballast at the pre-vacuum pump. |
| | Microscope does not vent. | No nitrogen. | Check nitrogen supply. |
| | | No compressed air. | Check compressed air supply. |
| Vac ready = OK is displayed abnormally fast. | Penning gauge has not been identified correctly. | Restart the microscope. | |

| Keyword | Symptom | Cause | Remedy |
|---------|---|---|--|
| | Gun vacuum is worse than $8 \text{ to } 9 \times 10^{-9}$ mbar. | The pumping capacity of the ion getter pump decreases in the course of time, thus deteriorating the gun vacuum. | If this does not solve the problem, contact your local ZEISS service representative. Refer to <i>Baking out the Gun Head</i> [▶ 102]. |

Optional Components

| Keyword | Symptom | Cause | Remedy |
|-----------|--|---|--|
| Flood Gun | No communication can be established with Smart-SEM | Loose cable connection. Flood Gun controller has an old firmware version. Flood Gun is not activated in FIB Configurator. | Check the cable connection. If this does not help, contact the ZEISS service representative. |

8.2 Overall System

8.2.1 Checking the CAN Communication

Purpose Checking the CAN Communication is useful if the microscope does not react to your commands anymore.

- In the Panel Configuration Bar, double-click CAN Communication.
→ The CAN Communication window is displayed.
- If any of the values is indicated as Yes, make sure that all cable connections between workstation and PC are plugged in correctly.
- If this does not help, reset the workstation as described in the instruction manual of the microscope.
INFO: If the problem persists, contact your ZEISS service representative.

8.3 Chamber

8.3.1 Initializing the Stage

Purpose If a stored stage position cannot be approached or if the stage does not move or does not move accurately, the stage needs to be initialized.

- Prerequisite**
- ✓ The specimen chamber has been evacuated, see *Loading the Specimen Chamber* [▶ 70]
 - ✓ Requires the Stage Initialise privilege
 - ✓ If there are any large specimens inside the chamber, remove them before initializing
- From the Menu Bar, select Stage > Stage Initialise.
→ The **Initialise Stage** window is displayed.
 - Confirm via Yes.
→ The stage initialization process takes a few minutes.

- INFO: If initialization of the stage does not solve the stage problem, contact your local ZEISS service representative.

8.3.2 Defining the Post Initialization Position of the Stage

Purpose You can configure the position to which the stage drives after the initialization procedure. Otherwise, the stage drives to the center position.

Prerequisite ✓ Requires the Supervisor privilege.

- From the Windows start menu, select SmartSEM > SmartSEM Administrator.
 - The SmartSEM Administrator Log on window is displayed.
- Enter user name and password.
- To confirm, click OK.
 - The SmartSEM Administrator window is displayed showing the user list.
- Click Column/Stage.
- In the Stage Post Initialisation Position input fields, enter the desired position. Alternatively, use the dual joystick to navigate to the desired position and click Set to current position.
- To activate the function, activate the Post Init. Posn Valid checkbox.

8.3.3 Changing the Joystick TV Angle


Purpose In TV mode (chamberscope), it can occur that dual joystick and stage seem to move to opposite directions. This is because the selected CCD camera is installed at a certain angle relative to the stage. Thus, the camera shows a side-inverted view. To remedy this, you need to change the joystick TV angle setting in the software.

Info

If you are working with two CCD cameras: The joystick TV angle can only be set for one CCD camera. When selecting the other CCD camera, you have to change the setting.

Prerequisite ✓ Requires the Supervisor privilege.

- From the Windows start menu, select SmartSEM > SmartSEM Administrator.
 - The SmartSEM Administrator Log on window is displayed.
- Enter user name and password.
- To confirm, click OK.
- The SmartSEM Administrator window is displayed showing the user list.
- Click Column/Stage.
- In the Stage Options section, double-click the Joystick TV Angle input field.
- Enter an angle depending on the installation location of the CCD camera.
 - If the CCD camera is installed at the back, enter 180°.
 - If the CCD camera is installed at the front, enter 0°.
 - If the CCD camera is installed at the side, enter 90°.



8.3.4 Replacing the Chamber Door Seal

Purpose Possible reasons for replacing the chamber door seal are the following:

- Chamber door does not close tightly
- Bad chamber vacuum

Overview The procedure contains the following steps:

- *Venting the Specimen Chamber* [▶ 70]
- *Replacing the O-ring* [▶ 101]
- *Evacuating the Specimen Chamber* [▶ 72]

8.3.4.1 Replacing the O-ring

WARNING

Suffocation hazard: Lack of oxygen

Gaseous dry nitrogen is used to vent the specimen chamber during specimen exchange. Inhaling nitrogen may cause unconsciousness.

- ▶ During specimen exchange, keep the chamber door open as short as possible.
- ▶ Do not inhale the air from within the specimen chamber.
- ▶ Ensure that the area around the microscope is sufficiently vented.
- ▶ If you begin to experience symptoms of asphyxia (for example: rapid breathing, loss of mental alertness and/or muscular coordination, depression of sensations, emotional instability, fatigue) leave the room immediately and inform the facility's safety officer.

NOTICE

Risk of property damage: Contamination caused by fingerprints

Contamination caused by fingerprints can lead to vacuum deterioration or prolonged pumping times.

- ▶ Always wear lint-free gloves when touching the specimen, specimen holder, or stage.

1. Carefully open the chamber door.
2. On the inside of the chamber door, remove the chamber door O-ring.
 - NOTICE Risk of property damage: If you use a metal tool to remove the O-ring, then you may damage the sealing surface. If necessary, then only use a plastic or wooden tool to remove the O-ring.**
3. Inspect the groove that holds the O-ring and remove any contamination.
4. Insert the new chamber door O-ring.
5. Close the chamber door.

8.3.5 Resetting the Touch Alarm

Purpose To prevent damage, a touch alarm is integrated in the microscope. If the specimen or the specimen holder touches the chamber walls, the detectors, or the objective lens, the stage is stopped immediately. An audible warning sounds and an on-screen message is displayed.

Prerequisite ✓ The EM Server shows the message WARNING Stage Touching.

1. To accept the warning, click OK.
2. Move the stage in the reverse direction away from the touch.

8.3.6 Checking the Water Flow and Temperature

1. In the Panel Configuration Bar, double-click Water Flow/Temperature.
 - The Water Flow/Temperature panel is displayed.
2. Check the entries.
 - If a value is critical, it is displayed in red.

8.4 Column

8.4.1 Baking out the Gun Head

Purpose The pumping capacity of the ion getter pump decreases in the course of time, thus deteriorating the gun vacuum. This can be remedied by an ion getter pump bakeout as a regular maintenance procedure.

Overview The procedure contains the following steps:

- *Switching off the Gun* [▶ 102]
- *Starting the Bakeout* [▶ 103]
- *Switching on the Gun* [▶ 73]

8.4.1.1 Switching off the Gun

Safety Information

NOTICE

Risk of property damage: Schottky field emitter

If the Schottky field emitter is switched on and off too frequently or inappropriately, its lifetime is reduced.

- ▶ Avoid switching off the gun during the working week.
- ▶ Use Standby mode for the weekend or a break of up to a week.
- ▶ When using the Standby mode, activate the Partial Vent on Standby function.

1. In the right part of the Status Bar, click **Gun: ✓** or **All: ✓**.
 - The pop-up menu for Vacuum, Gun and EHT activation is displayed.
2. Click Shutdown Gun.
3. Wait until the gun has ramped down.
 - This may take up to 5 minutes.

8.4.1.2 Starting the Bakeout

Info

You cannot work with the microscope while the bakeout procedure runs.

Safety Information**NOTICE****Risk of property damage: Hot surfaces during bakeout**

Parts of the enclosure in the upper range of the column may become hot during bakeout, particularly after a long bakeout cycle.

- ▶ Do not place any objects on the grids of the electron column during bakeout.
- ▶ After the bakeout procedure, let surfaces cool down before working around the column.
- ▶ Only advanced operators are allowed to perform the bakeout procedure.

- Prerequisite**
- ✓ Requires the Supervisor privilege and the user level Service
 - ✓ Only advanced operators are allowed to perform the bakeout procedure
1. In the **Panel Configuration Bar**, double-click **Bakeout**.
→ The **Bakeout** dialog is displayed.
 2. If the **Full service bakeout** checkbox is available, deactivate the **Full service bakeout** checkbox.
INFO: Full service bakeout includes column heating that may lead to column misalignment.
 3. From the **Bakeout** drop-down list, select a bakeout cycle.
For 2 hours heating / 1.5 hours cooling, select **Quick**.
For 8 hours heating / 1.5 hours cooling, select **Overnight**.
For 43 hours heating / 7 hours cooling, select **Weekend**.
For a cycle defined by the operator, select **User**.
 4. To start the bakeout procedure, click **Bakeout Start**.

8.4.2 Calibrating the Probe Current

Purpose This function enables you to automatically calibrate the probe currents within a few minutes.

Calibrating the probe current can be necessary in the following cases:

- Before performing analytical applications (e.g. EDX, WDX)
- After changing the extractor voltage
- To improve the accuracy of the set probe current values

Parts and Tools**Designation****Part no.**

Faraday cup

348342-8055-000

- Prerequisite**
- ✓ The microscope has a Gemini II column
1. Load the Faraday cup into the specimen chamber.
 2. Pump the specimen chamber.
 3. Switch on the electron beam.
 4. Set a magnification that allows transmission of the complete electron beam into the cavity of the Faraday cup through the aperture orifice.
 5. In the Panel Configuration Bar, double-click Probe Current Calibration.
→ The Probe Current Calibration window is displayed.
 6. Activate the Spot checkbox.

7. Click Cal I Probe.
8. To confirm, click Yes.
9. To store the calibration, click Save.
10. Deactivate the Spot checkbox.

8.5 Power Circuit

8.5.1 Checking the Position of the Circuit Breakers

Safety Information

NOTICE


Risk of property damage: Persisting electrical problems

Tripped circuit breakers may be a hint for an electrical problem in the microscope.

- ▶ If a circuit breaker keeps tripping, de-energize the microscope completely and contact your ZEISS service representative for assistance.

| No. | Value | Circuit |
|-----|-------|--------------------------|
| F1 | 10 A | Power supply unit |
| F2 | 10 A | PC, WDX, EDX, AUX 1–4 |
| F3 | 10 A | Airlock, pre-vacuum pump |
| F4 | 10 A | Internal heaters |

1. Check if one of the circuit breakers on the rear side of the plinth is tripped.



2. If one of the circuit breakers is tripped, push it upwards.

8.6 Detectors

8.6.1 Lubricating the Rod

The rod from the aSTEM and BSD4 detector mechanics needs to be lubricated once a year with TEM oil 300.

| Parts and Tools | Designation | Part no. |
|-----------------|------------------|-----------------|
| | TEM oil 300 | 000000-0484-955 |
| | Isopropanol | – |
| | Lint-free cloth | – |
| | Lint-free gloves | – |

CAUTION

Risk of injury: TEM oil 300

TEM oil 300 may be irritating to skin and eyes.

- ▶ Avoid contact with skin.
- ▶ Wear suitable gloves.
- ▶ In case of contact with the eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ After contact with the skin, wash immediately with plenty of water and soap.

CAUTION

Risk of injury: Isopropanol

Isopropanol is highly flammable and irritating to the eyes.

Vapors may cause drowsiness and dizziness.

- ▶ Wear suitable gloves.
- ▶ Keep away from sources of ignition.
- ▶ Do not smoke.
- ▶ In case of contact with the eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ Avoid contact with skin.
- ▶ Do not breathe vapor.

NOTICE

Risk of property damage: Unsuitable lubricants

When using unsuitable lubricants, the vacuum system may be contaminated.

- ▶ Only use TEM oil 300 for lubricating.

1. Retract the respective detector.
2. Clean the rod with isopropanol with a clean, lint-free cloth.
3. Spread some drops of TEM oil 300 across the rod.
Use a clean, lint-free cloth.

8.7 Focused Ion Beam

8.7.1 Checking the Lifetime of the Ion Source

Purpose The ion source is a consumable that is used up during operation. In order to keep track of the ion supply, the ion source lifetime needs to be checked regularly.

Info

When the lifetime of the ion source approaches 3000 μAh , contact the ZEISS service representative to have the ion source replaced.

1. In the **FIB Toolbar**, from the **FIB** drop-down list, select **FIB**.
 - The **FIB** tab of the **FIB Control Panel** is displayed.
 - The lifetime of the ion source is indicated as μAh .
2. Check the lifetime of the ion source.

8.7.2 Regenerating the Ion Source by Heating

Purpose From time to time, the gallium emitter has to be regenerated by heating. The heating procedure is used in following cases:

- If the suppressor voltage has reached its limits (i.e. 0 V or 2 kV) while a probe current of 2 μA cannot be maintained anymore
- If the ion source does not start emitting
- If the source is switched into Standby mode because the emission current exceeds the threshold defined in the SmartSEM Administrator in the Go to Standby if emission exceeds input field.

Info

Source heating can also be performed via Source Management. Source Management includes automatically triggered procedures to ensure an optimum performance of the ion source. These procedures are configured in the SmartSEM Administrator via FIB > FIB Source Management.

Safety Information

NOTICE

Risk of property damage: Ion source

Danger of damaging the ion source due to overheating. Excessive filament heating may reduce source life.

- ▶ Only heat as often as necessary.

Prerequisite ✓ Requires the Supervisor privilege.

1. In the Panel Configuration Bar, double-click FIB Ga Source Heating.
 - The FIB Ga Source Heating window is displayed.
2. Click Multi-Heat.
 - An automatic heating routine is performed.
3. If the ion source still does not start emitting, find a new extractor value. Refer to *Manually Finding a New Extractor Value* [▶ 107].

8.7.3 Manually Finding a New Extractor Value

Purpose If emission still fails after performing the multi-heating routine (see *Regenerating the Ion Source by Heating* [▶ 106]), you have to find a new extractor value.

In a first step, you can use the automatic FIB extractor adjustment functions:

- You can use Source Management. Source Management includes automatically triggered procedures to ensure an optimum performance of the ion source. These procedures are configured in the **SmartSEM Administrator** via **FIB > FIB Source Management**.
- You can trigger the extractor adjustment via **FIB Control Panel > FIB > Extr. Adjust**.

If the automatic FIB extractor adjustment functions fail, you have to find a new extractor value manually.

1. In the **FIB Toolbar**, from the **FIB** drop-down list, select **FIB**.
→ The **FIB** tab of the **FIB Control Panel** is displayed.
2. Click **Standby**.
→ The FIB gun is set to **Standby** mode.
3. To disable the suppressor regulation, deactivate the **Regulate** checkbox.
4. Set the **FIB Suppressor Target** to 1200 V.
5. Select the **Options** tab.
6. Activate the **FIB Extractor Manual Override** checkbox.
7. Move the slider **FIB Manual Extractor Target** slowly until FIB emission starts.
8. Adjust the **FIB Manual Extractor Target** until emission is approximately 2 μA .
9. Leave the source in this condition for about 5 minutes, adjust the extractor voltage if necessary to keep the emission at 2 μA .
10. Click **Save Value**.
11. Deactivate the **FIB Extractor Manual Override** checkbox.
12. Select the **FIB** tab.
13. To enable the suppressor regulation, activate the **Regulate** checkbox.

Info

During the process of finding a new extractor target, the FIB gun will be automatically switched to the **On** state.

9 Shutdown and Disposal

9.1 Putting the Microscope out of Operation

If the microscope will not be used for an extended period of time, e.g. several months, it should be put out of operation.

Contact your local ZEISS service representative to have the microscope put out of operation.

9.2 Transport and Storage

WARNING

Risk of injury: Tilting hazard when removing the microscope from the crate

When removing the microscope from the wooden crate, it can tilt and crush a person.

- ▶ Use a forklift to remove the microscope from the wooden crate.

CAUTION

Risk of injury: Crushing hazard when lowering the microscope

The microscope and its components are heavy. When the load is lowered during transport and positioning, body parts can be crushed.

- ▶ Maintain a safe distance.
- ▶ Do not walk or place your hands or feet under the load while it is being lowered.
- ▶ Wear safety shoes and gloves.

CAUTION

Risk of injury or property damage: GIS precursors

If a gas injection system (GIS) is used, the precursors may be explosive, reactive, toxic or irritant when they come into contact with the environment or a person.

- ▶ Do not remove or install a GIS reservoir. Contact your local ZEISS service representative to have a reservoir removed or installed.
- ▶ The microscope may only be shipped if all precursors have been removed.
- ▶ Precursors have to be shipped separately in a special transport container. For shipment always enclose a print-out of the respective material safety data sheet (MSDS).
- ▶ Store and transport precursors only within the prescribed temperature ranges.
- ▶ For further information, refer to the GIS instruction manual.

NOTICE**Risk of property damage: Damage during transport**

Sensitive components of the microscope can get damaged during transport.

- ▶ The microscope may only be transported in air-suspended vehicles.
- ▶ Moving parts must be secured during transport to prevent them from slipping or tipping over.
- ▶ Install shock/tilt watches.
- ▶ Avoid rocking the crates back and forth.
- ▶ Devices for transporting the microscope must be rated to handle its full weight and dimensions. Note the weight information on the package and on the shipping document.
- ▶ Check that none of the items has been damaged during shipment.
- ▶ Otherwise contact your local ZEISS service representative.

Packaging of the Microscope

The microscope is delivered in two wooden crates:

- Microscope plinth

Wrapped with recyclable polyethylene-foil and shipped in a reusable wooden crate.

Dimensions and weight of crate:

1550 × 1550 × 2200 mm³ (W × D × H), appr. 1300 kg

- Microscope console and accessories

Console, valve, damper, monitors, cables, pipes etc. are wrapped with recyclable polyethylene-foil or packed in separate cartons and shipped in a reusable wooden crate.

Dimensions and weight of crate:

1720 × 1100 × 980 mm³ (W × D × H), appr. 300 kg

Guidelines for Unpacking the Microscope

Due to the heavy weight of microscopes, a forklift has to be used to remove the microscopes from the wooden crate:

- The forklift used must have a sufficient load capacity.
- Refer to the weights of the microscope stated in this chapter.

Guidelines for Transporting the Unpacked Microscope

A hand pallet truck has to be used to move the microscope:

- Ensure all hallways and corners are wide enough to be passed with the hand pallet truck.

Allowable Conditions during Storage and Transport

The packed microscope has to be stored in a dry place.

Allowable temperature during storage and transport:

Microscope: between $-10\text{ }^{\circ}\text{C}$ and $+70\text{ }^{\circ}\text{C}$

If a GIS system is part of the microscope, the GIS precursors have to be transported and stored within the prescribed temperature ranges.

Allowable temperature during storage and transport:

| | | |
|------------|------------------|---|
| Precursor: | Gold | $-30\text{ }^{\circ}\text{C}$ (at least $< 8\text{ }^{\circ}\text{C}$) |
| | Platinum | $< 5\text{ }^{\circ}\text{C}$ |
| | XeF ₂ | $< 5\text{ }^{\circ}\text{C}$ |

All other reservoirs are allowed to be stored and transported without cooling. However, ZEISS recommends to cool all reservoirs during transport.

9.3 Disposal**9.3.1 Disposing of Solid Waste (Consumables)**

The operator must ensure that solid waste (consumables) is disposed of and recycled in a responsible manner.

| Description | Material | Disposal |
|------------------------|---|---|
| Schottky field emitter | Tungsten, ceramics | Used up emitters may be returned to the emitter manufacturer. |
| Apertures | Platinum, iridium, gold | Very small amounts. May be disposed of in accordance with local/regional regulations. |
| Ion source (gallium) | Gallium | To be returned to ZEISS. |
| GIS Precursors | For information regarding the disposal of precursors, refer to the instruction manual of the GIS. | |

9.3.2 Disposing of the Microscope

The operator must ensure that waste products are disposed of and recycled in a responsible manner.

Refer to EC directive 2012/19/EC on waste electrical and electronic equipment (WEEE).

The microscope consists of several modules. Be careful to separate the materials properly when you dispose of the microscope.

- Materials: e.g. metals, non-metals, composite materials, process materials
- Electronic scrap material: e.g. transformers, circuit boards, cables

Comply with national and regional waste disposal regulations.

10 Technical Data and Conformity

10.1 Product Specification

Electron Optics

| Parameter | Description |
|-----------------------------|--|
| Resolution | GEMINI II column in High Resolution configuration (max. probe current 40 nA) |
| | Optimum working distance Without using Tandem Decel: |
| | <ul style="list-style-type: none"> ▪ 0.9 nm at 15 kV ▪ 1.6 nm at 1 kV |
| | Including Tandem Decel option: |
| | GEMINI II column in High Current configuration (max. probe current 100 nA) |
| | Optimum working distance Without using Tandem Decel: |
| | <ul style="list-style-type: none"> ▪ 0.9 nm at 15 kV ▪ 1.8 nm at 1 kV |
| | Including Tandem Decel option: |
| Acceleration voltage | Range: 20 V to 30 kV |
| | Adjustment: Continuously variable in 10 Volt steps |
| Probe current | High Resolution configuration: 10 pA to 40 nA High Current configuration: 10 pA to 100 nA |
| | Probe current regulation via a double condenser system |
| Magnification | Range: 12x to 2,000,000x referenced to Polaroid image format |
| Electron source | Filament: Schottky field emitter |
| | Automatic emitter run-up: Safe controlled run-up to the target emitter conditions |
| Objective lens | Type: Patented GEMINI II electromagnetic/electrostatic objective lens system (68° conical final lens) with water cooling for best thermal stability and reproducibility |
| Focus | Working distance: Range from 1 to 50 mm, depending on acceleration voltage |
| | Focus compensation: Automatic compensation to minimize focus changes over the entire acceleration voltage and probe current range |
| | Dynamic focus: For correction of focus on tilted specimens |
| | Focus wobble: For assistance in aperture alignment, with adjustable amplitude and speed |

| Parameter | Description |
|--------------------------------|--|
| | Rotation compensation: Automatic correction of apparent image rotation with changes in working distance |
| Beam shift | For precise adjustment of image position at high magnifications Width: 200 µm (±100 µm) at 20 kV and WD = 8.5 mm |
| Tandem decel (optional) | Tandem decel module to apply negative bias to the specimen. Bias voltage can be set to 1 kV, 2 kV, 3 kV, 5 kV and continuously varied between 50 V and 100 V in steps of 1 V. Airlock enabled specimen holder designed to ensure optimum electrical field geometry included Bias wobble for assistance with precise tilt alignment of specimen |

Focused Ion Beam Column

| Parameter | Description |
|-----------------------------|---|
| Ion Source | Type: UHV, with gallium liquid metal ion source (Ga-LMIS) Source Life: 3000 hours at 1 µA of emission current Isolation: Automatically controlled valve for source isolation |
| FIB resolution | 3 nm at 30 kV |
| Probe current | 1 pA to 100 nA |
| Magnification | 300x to 500,000x |
| Acceleration Voltage | Range: 0.5 to 30 kV Adjustment: Continuously variable in 10 Volt steps |
| Imaging/Patterning | Maximum field of view: 580 × 580 µm User beam shift: ±15 µm |
| Dwell time | Minimum dwell time: 25 ns Maximum dwell time: 1 s |
| Lenses | Type: Two electrostatic lenses |

Specimen Chamber and Stage

| Parameter | Description | |
|------------------------------------|--|---|
| | Crossbeam 550 | Crossbeam 550 L |
| Specimen chamber dimensions | <ul style="list-style-type: none"> ▪ 330 mm inner diameter ▪ 270 mm height | <ul style="list-style-type: none"> ▪ 520 mm inner diameter ▪ 307 mm height |
| Free accessory ports | 18 for EDS, EBSD, SIMS, manipulators etc. | 22 free ports for optional accessories such as load lock, STEM, 4QBSD, GIS, EDS, WDS, EBSD, CL, flood gun, cryo transfer, SIMS, manipulator, etc. |

| Parameter | Description | | | | | | | | | | | |
|------------------------------------|---|-----------------|----------------|----------------|-------------|-------------|-------------|-------------|---------------------|----------------------|------------------------------|------------------------------|
| | Crossbeam 550 | Crossbeam 550 L | | | | | | | | | | |
| Analytical working distance | 5 to 8.5 mm | | | | | | | | | | | |
| Coincidence point | 5 mm (SEM), 12 mm (FIB) | | | | | | | | | | | |
| Column arrangement | The angle between electron-optical and focused ion beam-column is 54°. | | | | | | | | | | | |
| Specimen stage | <p>Type: 6-axes motorized super-eucentric, controlled via the SmartSEM user interface, operated by a dual joystick control box</p> <p>Mounting: Drawer-type door</p> <p>Movements:</p> <table border="0"> <tr> <td>▪ X/Y = 100 mm</td> <td>▪ X/Y = 153 mm</td> </tr> <tr> <td>▪ Z = 50 mm</td> <td>▪ Z = 50 mm</td> </tr> <tr> <td>▪ M = 13 mm</td> <td>▪ M = 20 mm</td> </tr> <tr> <td>▪ Tilt = -4° to 70°</td> <td>▪ Tilt = -15° to 70°</td> </tr> <tr> <td>▪ Rotation = 360° continuous</td> <td>▪ Rotation = 360° continuous</td> </tr> </table> <p>INFO: The movements may be reduced by specimen size, operating conditions, and accessories attached.</p> <p>Accessory ports: Two accessory ports on the stage door are provided</p> <p>Specimen weight: Up to 0.5 kg</p> <p>Specimen current monitor with integrated touch alarm (audible touch alarm warning with on-screen message)</p> <p>Specimen mounts: One carousel 9x specimen holder for 13 mm diameter stubs included in base tool configuration; various specimen holders available as option</p> | | ▪ X/Y = 100 mm | ▪ X/Y = 153 mm | ▪ Z = 50 mm | ▪ Z = 50 mm | ▪ M = 13 mm | ▪ M = 20 mm | ▪ Tilt = -4° to 70° | ▪ Tilt = -15° to 70° | ▪ Rotation = 360° continuous | ▪ Rotation = 360° continuous |
| ▪ X/Y = 100 mm | ▪ X/Y = 153 mm | | | | | | | | | | | |
| ▪ Z = 50 mm | ▪ Z = 50 mm | | | | | | | | | | | |
| ▪ M = 13 mm | ▪ M = 20 mm | | | | | | | | | | | |
| ▪ Tilt = -4° to 70° | ▪ Tilt = -15° to 70° | | | | | | | | | | | |
| ▪ Rotation = 360° continuous | ▪ Rotation = 360° continuous | | | | | | | | | | | |
| Plasma cleaner (optional) | <p>Integrated plasma cleaner:</p> <p>Chamber mounted plasma cleaner for removal of hydrocarbon contamination from both specimen and chamber; integrated software control for user defined cleaning cycles without user interaction after starting the cleaning process</p> | | | | | | | | | | | |

Detection System

| Parameter | Description |
|-------------------------|--|
| Detection system | <p>Detection system equipped for parallel detection and signal processing of multiple detector channels:</p> <ul style="list-style-type: none"> ▪ Up to 12 detector channels and up to 3 TV inputs ▪ Parallel detection, processing, and display of up to four channels possible ▪ 2 signal detector inputs can be mixed for enhanced image information |
| InLens detectors | InLens SE detector: |

| Parameter | Description |
|---------------------------------|---|
| | High efficiency annular scintillator detector mounted in GEMINI column with optically coupled photomultiplier. |
| | InLens EsB detector (optional): Column-mounted high efficiency scintillator detector with optically coupled photomultiplier for detection of energy and angle selective backscattered electrons. Filtering grid adjustable from 0 V to –1.5 kV for contrast adjustment. |
| Chamber detectors | SE detector: Everhart-Thornley SE detector with optically coupled photomultiplier; collector bias adjustable from –250 to +400 V |
| | SESI detector (optional): Combined Secondary Electron Secondary Ion (SESI) detector based on a scintillator photomultiplier system; easy change between secondary ion and secondary electron mode by converting the electrode voltage; replaces the SE detector |
| | CL detector (optional): Cathodoluminescence (CL) chamber detector |
| | aSTEM detector (optional): Pneumatically retractable multi-mode annular Scanning Transmission Electron Microscopy (aSTEM) detector; enables bright field (BF), dark field (DF), and high angle annular dark field (HAADF) detection |
| | aBSD/BSD detector (optional): Pneumatically retractable 5 or 6 segment multi-mode solid state BSE detector; enables materials contrast, crystal orientation, and topographic imaging |
| Chamber camera | Color CCD camera with white-light illumination and IR illumination |
| Specimen current monitor | 6 range auto ranging for precise current measurement in the range of 1 pA to 10 μ A |

For more details, refer to the document Product Specification.

10.2 Installation Requirements

For a complete list of the installation requirements, refer to the document Installation Requirements.

Location Requirements

| Parameter | Requirement |
|------------------------------|------------------------------|
| Installation site | Exclusively inside buildings |
| Recommended room size | Min. 4.0 m × 5.0 m × 2.3 m |
| Service area | Min. 1.0 m at each side |
| Entrance | Min. 1.1 m wide |
| Hallways | Min. 1.3 m wide |

| Parameter | Requirement |
|-----------------------|------------------------------|
| Corners | Min. 1.5 m wide |
| Transport ways | Free of staircases |
| Installation category | II |
| Floor stability | > 1000 kg/m ² |
| Parameter | Requirement |
| Installation site | Exclusively inside buildings |
| Recommended room size | Min. 3.5 m × 5.0 m × 2.3 m |
| Service area | Min. 1.0 m at each side |
| Entrance | Min. 0.8 m wide |
| Hallways | Min. 1.0 m wide |
| Corners | Min. 1.2 m wide |
| Transport ways | Free of staircases |
| Installation category | II |
| Floor stability | > 1000 kg/m ² |
| Parameter | Requirement |
| Installation site | Exclusively inside buildings |
| Recommended room size | Min. 3.6 m × 5.0 m × 2.3 m |
| Service area | Min. 0.8 m at each side |
| Entrance | Min. 0.8 m wide |
| Hallways | Min. 1.0 m wide |
| Corners | Min. 1.2 m wide |
| Transport ways | Free of staircases |
| Installation category | II |
| Floor stability | > 1000 kg/m ² |
| Parameter | Requirement |
| Installation site | Exclusively inside buildings |
| Recommended room size | Min. 3.6 m × 4.0 m × 2.3 m |
| Service area | Min. 0.8 m at each side |

| Parameter | Requirement |
|-----------------------|--------------------------|
| Entrance | Min. 0.8 m wide |
| Hallways | Min. 1.0 m wide |
| Corners | Min. 1.2 m wide |
| Transport ways | Free of staircases |
| Installation category | II |
| Floor stability | > 1000 kg/m ² |

Exhaust Line

An exhaust line is required to remove the waste gas of the pre-vacuum pump and to transmit it to the outside.


If toxic chemicals or biological specimens are used an exhaust line is recommended to remove the waste gas of the pre-vacuum pump and to transmit it to the outside.

Electrical Supplies

| Parameter | Requirement |
|---|--|
| Nominal AC voltage | 208–230 VAC ($\pm 10\%$), L1/N(L2)/PE The provided electrical connection must be in accordance with the applicable electrical codes for the country of installation. In order to avoid disturbance from other installed machines, ZEISS recommends using a separate power connection to the main distribution panel. |
| Protection class | Class I |
| Nominal frequency | 50–60 Hz |
| Momentary interruption | Less than a half cycle |
| System connection | The microscope is delivered with a power cord 3 × AWG10 UL-style (3 m long) and a CEE-connector (Type 32 A-6 h; 200 to 230 V, 2P+PE). Alternatively the microscope can be installed directly to a switchable power distribution terminal which can be secured against accidental re-activation. |
| Power consumption | Max. 3.7 kVA, dependent on accessories |
| Maximum current | 16 A |
| Circuit breaker (at house installation) | 25 A (type K) |
| Ampere interrupting capacity (AIC) | Min. 10,000 A rms |
| Protective ground | High leakage currents are present in the microscope. Therefore, the microscope has to be connected to an equipotential bonding bar. An exclusive grounding connection to earth must be provided, i.e. the grounding terminal must not be common to other electrical equipment. A grounding wire AWG10 is delivered with the microscope. Cross section: > 4 mm ² |

| Parameter | Requirement |
|--|---|
| | Ground resistance: < 0.1 Ω |
| Parameter | Requirement |
| Nominal AC voltage | 230 V _{AC} 1/N/PE or 208 V _{AC} 2/PE The provided electrical connection must be in accordance with the applicable electrical codes for the country of installation. In order to avoid disturbance from other installed machines, ZEISS recommends using a separate power connection to the main distribution panel. |
| Protection class | Class I |
| Nominal frequency | 50–60 Hz |
| Momentary interruption | Less than a half cycle |
| System connection | The microscope is delivered with a 3 m long power cord that is equipped with a 200–250 V _{AC} CEE MALE PLUG 2P3W 6h 16A (blue) according to IEC 60309. The supplied connection cable must not be replaced with another one. Otherwise, the conformity with listed standards becomes invalid. The building installation should provide the corresponding 200–250 V _{AC} CEE FEMALE RECEPTACLE 2P3W 6h 16A (blue) with the correct wiring 1/N/PE or 2/PE and the desired approvals of the country used. To avoid disturbance from other installed machines Carl Zeiss Microscopy recommends to use a separate power connection to the main distribution panel. |
| Power consumption | Max. 3 kVA, dependent on accessories |
| Circuit breaker (at house installation) | 16 A Type C |
| Ampere interrupting capacity (AIC) | Min. 10,000 A _{rms} |
| Protective ground | High leakage currents are present in the microscope. Therefore, the microscope has to be connected to a separate protective ground. An exclusive grounding connection to earth must be provided as part of the building installation, i.e. a grounding screw terminal (\varnothing 8 mm) which is directly connected to the PE of the FEMALE RECEPTACLE as short as possible. (see picture) This grounding connection must not be common to other electrical equipment. |


| Parameter | Requirement |
|-----------|-------------|
|-----------|-------------|



A grounding wire AWG10 (≥ 5 m) is delivered with the microscope. It serves as connection between the grounding screw terminal (Ø 8 mm) and the microscope (Ø 6 mm).

Cross section: Min. AWG10 (between grounding screw terminal and PE of the CEE FEMALE RECEPTACLE)²

| Parameter | Requirement |
|--|--|
| Nominal AC voltage | 120 V _{AC} (single phase) or 230 V _{AC} (single phase) The provided electrical connection must be in accordance with the applicable electrical codes for the country of installation. In order to avoid disturbance from other installed machines, ZEISS recommends using a separate power connection to the main distribution panel. Details see System Connection |
| Protection class | Class I |
| Nominal frequency | 50–60 Hz |
| Momentary interruption | Less than a half cycle |
| System connection | The microscope is delivered with a power cord 3 m long equipped with a CEE plug 1/N/PE according to IEC 60309 depending on the customer's site supply: either for 120 V _{AC} a CEE MALE PLUG 2P3W 4h 20A (yellow) or for 230 V _{AC} a CEE MALE PLUG 2P3W 6h 16A (blue). The building installation should provide the corresponding CEE socket 1/N/PE with desired approvals of the country used: either for 120 V _{AC} a CEE FEMALE RECEPTACLE 2P3W 4h 20A (yellow) or for 230 V _{AC} a CEE FEMALE RECEPTACLE 2P3W 6h 16A (blue). To avoid disturbance from other installed machines Carl Zeiss Microscopy recommends to use a separate power connection to the main distribution panel. |
| Power consumption | Max. 3 kVA, dependent on accessories |
| Circuit breaker (at house installation) | 16 A Type C |
| Ampere interrupting capacity (AIC) | Min. 10,000 A _{rms} |

| Parameter | Requirement |
|--------------------------|---|
| Protective ground | <p>High leakage currents are present in the microscope. Therefore, the microscope has to be connected to a separate protective ground.</p> <p>An exclusive grounding connection to earth must be provided as part of the building installation, i.e. a grounding screw terminal (Ø 8 mm) which is directly connected to the PE of the FEMALE RECEPTACLE as short as possible.</p> <p>This grounding connection must not be common to other electrical equipment.</p>  <p>A grounding wire AWG10 (≥ 5 m) is delivered with the microscope. It serves as connection between the grounding screw terminal (Ø 8 mm) and the microscope (Ø 6 mm).</p> <p>Cross section: Min. AWG10 (between grounding screw terminal and PE of the CEE FEMALE RECEPTACLE)²</p> |

Cooling Water

Major components of the microscope such as electron-optic lenses, parts of the electronics and the turbo molecular pump are water-cooled. Any cooling solution has to fulfil the following requirements.

| Parameter | Requirement |
|--------------------------|---|
| Water flow rate | 60–70 l/h |
| Pressure | Adjustable up to 3 bar |
| Water temperature | 20–22 °C |
| Stability | 0.5 °C/10 min |
| Heat dissipation | 1 kW |
| Connection hose | 6 mm inside diameter. Two pieces 10 m each are delivered. |

| Parameter | Requirement |
|--------------------------|-------------|
| Water flow rate | > 1.5 l/min |
| Pressure | 2–3 bar |
| Water temperature | 18–22 °C |

| Parameter | Requirement |
|-----------------------|---|
| Stability | 0.5 °C/10 min |
| Heat dissipation | 1 kW |
| Connection hose | 8 mm inside diameter. One 25 m roll is delivered with the microscope. |
| Instrument connection | Quick exchange connectors. Two are delivered with the microscope. |

Nitrogen

Gaseous dry nitrogen is used to vent the specimen chamber during specimen exchange. The nitrogen can be taken either from a gas cylinder or from an in-house supply system.

The connection must be equipped with an appropriate pressure reducer and a shut-off valve that is secured against accidental re-activation.

| Parameter | Requirement |
|-----------------|---|
| Flow rate | Approx. 40 l/min for ventilation of specimen chamber with chamber door open |
| Pressure | 0.30–0.35 MPa (3.0–3.5 bar) |
| Quality | 4.6 with nitrogen content > 99.996 % |
| Connection hose | 6 mm inside diameter. 10 m are delivered with the microscope. |

| Parameter | Requirement |
|-----------------------|--|
| Flow rate | Approx. 3 l/min for ventilation of specimen chamber with chamber door open |
| Pressure | 0.2–3.3 bar |
| Quality | 4.6 with nitrogen content > 99.996 % |
| Connection hose | 4 mm inside diameter. 10 m are delivered with the microscope. |
| Instrument connection | Quick exchange connector. One is delivered with the microscope. |

Compressed Air

Compressed air is used to operate several valves and the auto leveling system.

The necessary compressed air can be either generated by a compressor (part no. 345596-0000-000) or taken from a gas cylinder or from an in-house supply system.

The connection must be equipped with an appropriate pressure reducer and a shut-off valve that can be secured against accidental re-activation.

| Parameter | Requirement |
|-------------------|---|
| Typical flow rate | Approx. 12 l/min at 0.6 MPa pressure during air leveling system inflation |
| Pressure | 0.6–0.8 MPa (6–8 bar) |
| Quality | Oil-free |
| Connection hose | 6 mm inside diameter. 10 m are delivered with the microscope. |

| Parameter | Requirement |
|-----------------------|--|
| Typical flow rate | Less than 1 l/min during normal operation |
| Pressure | 0.6–0.8 MPa (6–8 bar) |
| Quality | Oil-free |
| Connection hose | 4 mm inside diameter. 10 m of pipe are delivered with the microscope. |
| Instrument connection | Quick exchange connector. One is delivered with the microscope. INFO: Due to acoustic noise and vibrations the compressor – if used – should be installed in a separate room. |

Environmental Requirements

| Parameter | Requirement |
|---|--|
| Ambient temperature | Appr. 21 ± 4 °C |
| Stability of ambient temperature | 0.5 °C/h |
| For long-running experiments: Long-term stability of ambient temperature | 2 °C/24 h |
| Relative humidity | Less than 65 % |
| Altitude | Max. 2000 m above sea level to guarantee an undisturbed operation |
| Pollution degree | 2 |
| Electrical field | The microscope is a class A device (industrial). The microscope is designed to operate in a controlled electromagnetic environment. This means that devices with RF transmitters such as mobile phones or DECT phones must not be used in close proximity. |
| Vibrations | Up to 10 Hz: less than 5 $\mu\text{m/s}$ 10–60 Hz: less than 10 $\mu\text{m/s}$ Above 60 Hz: less than 14 $\mu\text{m/s}$ |
| Magnetic stray fields | Measured in time domain and at specimen chamber height (1.0–1.5 m). DC component: 0.5 mG / 5 min or less AC component: less than 1 mG peak to peak between 10 Hz and 1 kHz |
| Acoustic noise | Up to 120 Hz: less than 52 dB 120–450 Hz: less than 43 dB Above 450 Hz: less than 47 dB |


| Parameter | Requirement |
|---------------------|---------------------|
| Ambient temperature | Appr. 21 ± 4 °C |

| Parameter | Requirement |
|---|--|
| Stability of ambient temperature | 0.5 °C/h |
| For long-running experiments: Long-term stability of ambient temperature | 2 °C/24 h |
| Relative humidity | Less than 65 % |
| Altitude | Max. 2000 m above sea level to guarantee an undisturbed operation |
| Pollution degree | 2 |
| Electrical field | The microscope is a class A device (industrial). The microscope is designed to operate in a controlled electromagnetic environment. This means that devices with RF transmitters such as mobile phones or DECT phones must not be used in close proximity. |
| Vibrations | Horizontal vibrations (in x/y-direction) Up to 10 Hz: less than 5 µm/s 10–60 Hz: less than 10 µm/s Above 60 Hz: less than 14 µm/s Vertical vibrations (in z-direction) Up to 10 Hz: less than 4 µm/s 10–60 Hz: less than 14 µm/s Above 60 Hz: less than 20 µm/s |
| Magnetic stray fields | Measured in time domain and at specimen chamber height (1.0–1.5 m). DC component: 0.5 mG / 5 min or less AC component: less than 1 mG peak to peak between 10 Hz and 1 kHz |
| Acoustic noise | Up to 120 Hz: less than 52 dB 120–450 Hz: less than 43 dB Above 450 Hz: less than 47 dB |
| Parameter | Requirement |
| Ambient temperature | Appr. 21±4 °C |
| Stability of ambient temperature | 0.5 °C/h |
| For long-running experiments: Long-term stability of ambient temperature | 2 °C/24 h |
| Relative humidity | Less than 65 % |
| Altitude | Max. 2000 m above sea level to guarantee an undisturbed operation |


| Parameter | Requirement |
|---|--|
| Pollution degree | 2 |
| Electrical field | The microscope is a class A device (industrial). The microscope is designed to operate in a controlled electromagnetic environment. This means that devices with RF transmitters such as mobile phones or DECT phones must not be used in close proximity. |
| Vibrations | Horizontal vibrations (in x/y-direction) Up to 10 Hz: less than 0.04 mm/s 10–20 Hz: less than 0.17 mm/s 20–70 Hz: less than 0.30 mm/s Above 70 Hz: less than 20 mm/s Vertical vibrations (in z-direction) Up to 8 Hz: less than 0.03 mm/s 8–45 Hz: less than 0.15 mm/s Above 45 Hz: less than 2.00 mm/s |
| Magnetic stray fields | Less than 3 mG peak to peak between 10 Hz and 1 kHz (Sigma 300) Less than 1 mG peak to peak between 10 Hz and 1 kHz (Sigma 500) |
| Acoustic noise | Sigma 300: Less than 53 dB for frequencies up to 200 Hz Less than 42 dB for frequencies from 200 up to 300 Hz Less than 50 dB for frequencies higher than 300 Hz Sigma 500: Less than 50 dB for frequencies up to 200 Hz Less than 40 dB for frequencies from 200 up to 400 Hz Less than 45 dB for frequencies higher than 400 Hz |
| Parameter | Requirement |
| Ambient temperature | 20-30 °C It is recommended that for operator comfort the room temperature is maintained between 22 °C and 24 °C. |
| Stability of ambient temperature | 0.5 °C/h |
| For long-running experiments: Long-term stability of ambient temperature | 2 °C/24 h |
| Relative humidity | Less than 65 % |
| Altitude | Max. 2000 m above sea level to guarantee an undisturbed operation |
| Pollution degree | 2 |
| Electrical field | The microscope is a class A device (industrial). The microscope is designed to operate in a controlled electromagnetic environment. This means that devices with RF transmitters such as mobile phones or DECT phones must not be used in close proximity. |

| Parameter | Requirement |
|------------------------------|---|
| Vibrations | Less than 6 $\mu\text{m}/\text{sec}$ rms from 0–30 Hz Less than 12 $\mu\text{m}/\text{sec}$ rms above 30 Hz |
| Magnetic stray fields | Max. 3×10^{-7} T (peak to peak) = 3 mG (peak to peak) |
| Acoustic noise | Less than 53 dB for frequencies up to 200 Hz Less than 42 dB for frequencies from 200 up to 300 Hz Less than 50 dB for frequencies higher than 300 Hz |

10.2.1 Layout and Connections



| | |
|--|--|
| 1 Static vibration damper | A Mains power supply 208...230 V/25 A, L1/N(L2)/PE |
| 2 Pre-vacuum pump | B Equipotential bonding bar |
| 3 Computer workplace | C Pressure reducer (water, nitrogen, compressed air) |
| 4 Emergency Off (EMO) button (optional) | D Main shut-off valves |
| | E Water supply |
| | F Water runback |
| | G Nitrogen supply |
| | H Compressed air supply |
| | I Exhaust line |



1 Static vibration damper

2 Pre-vacuum pump

3 Plinth and column

4 Computer workplace

A Mains power supply 208...230 V,
max. 16 A


B Equipotential bonding bar, safety
earth

C Vacuum exhaust

D Chilled water supply

E Dry compressed air

F Dry nitrogen supply



1 Static vibration damper

A Mains power supply
120 V, 50–60 Hz, 16 A, 1/N/PE or
230 V, 50–60 Hz, 16 A, 1/N/PE

2 Pre-vacuum pump

B Equipotential bonding bar /
grounding screw terminal (Ø 8 mm)

3 Pre-vacuum pump 2 (optional)

C Pressure reducer (nitrogen)

4 Computer workplace

D Main shut-off valve

5 Emergency Off (EMO) box (optional)

E Nitrogen supply

6 Emergency Off (EMO) button (optional)

F Exhaust line (optional)

10.2.2 System Layout




Fig. 37: System layout

- 1** Plinth + column + Ion-sculptor FIB (+ 200 mm airlock)
- 2** Table + PC
- 3** Static damping block
- 4** Pre-vacuum pump
- 5** FIB rack

| Description | Size (mm) appr. | Distribution of load (kg) | Footprints |
|---|--|---------------------------|-------------|
| Plinth + column + Ion-sculptor FIB | 940 × 1070 × 1800 | 4 × 250 | 4 × Ø 80 mm |
| Plinth + column + Ion-sculptor FIB + 200 mm airlock | 1080 × 1070 × 1800 | 4 × 260 | 4 × Ø 80 mm |
| Table + PC | 1650 × 1000 × 838 or 1100 × 1000 × 838 | 4 × 24.3 | 4 × Ø 50 mm |

| Description | Size (mm) appr. | Distribution of load (kg) | Footprints |
|---|------------------|---------------------------|-----------------|
| Static damping block | 180 × 180 × 160 | 1 × 12.0 | 180 mm × 180 mm |
| Pre-vacuum pump | 432 × 265 × 295 | 1 × 24.5 | 200 mm × 180 mm |
| FIB rack | 560 × 810 × 1600 | 4 × 50 | on wheels |
| Chiller (optional, water- or air-cooled)* | 530 × 640 × 740 | 4 × 22.5 | on wheels |
| Compressor* | 281 × 139 × 239 | 1 × 25.0 | 281 mm × 139 mm |

*Actual measurements are country dependent (different manufacturers).

| Description | Size (mm) appr. | Distribution of load (kg) | Footprints |
|--|---|---------------------------|-----------------|
| Plinth + column + Ion-sculptor FIB | 940 × 1070 × 1800 (Crossbeam 550) | 4 × 250 | 4 × Ø 80 mm |
| | 940 × 1070 × 1880 (Crossbeam 550 L) | 4 × 297.5 | |
| Plinth + column + Ion-sculptor FIB + 200 mm airlock | 1080 × 1070 × 1800 (Crossbeam 550) | 4 × 260 | 4 × Ø 80 mm |
| | 1080 × 1070 × 1880 (Crossbeam 550 L) | 4 × 305 | |
| Table + PC | 1650 × 1000 × 838 | 4 × 24.3 | 4 × Ø 50 mm |
| | or 1100 × 1000 × 838 | | |
| Static damping block | 180 × 180 × 160 | 1 × 12.0 | 180 mm × 180 mm |
| Pre-vacuum pump | 432 × 265 × 295 | 1 × 24.5 | 200 mm × 180 mm |
| FIB rack | 560 × 810 × 1600 | 4 × 50 | on wheels |
| Chiller (optional, water- or air-cooled)* | 530 × 640 × 740 | 4 × 22.5 | on wheels |
| Compressor* | 281 × 139 × 239 | 1 × 25.0 | 281 mm × 139 mm |

*Actual measurements are country dependent (different manufacturers).

10.3 Declaration of Conformity

Denomination Focused Ion Beam - Scanning Electron Microscope (FIB-SEM)

Model Crossbeam 550/550L

Manufacturer Carl Zeiss Microscopy GmbH
Carl-Zeiss-Str. 22
73447 Oberkochen
Germany

We declare under our sole responsibility that the above machinery and the optional accessories fulfil all the relevant provisions of the following EC Directives:

- **2006/42/EC** – Machinery Directive
- **2014/30/EU** – Electromagnetic Compatibility

Applied harmonized standards:

- **EN ISO 12100:2010** – Safety of machinery – General principles for design – Risk assessment and risk reduction
- **EN ISO 13849-1:2008** – Safety of machinery – Safety related parts of control systems – Part 1: General principles for design
- **EN 60204-1:2006** – Safety of machinery – Electrical equipment of machines – Part 1: General requirements
- **EN 61010-1:2010** – Safety requirements for electrical equipment for measurement, control, and laboratory use
- **EN 61326-1:2013** – Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements

CE Marking The CE conformity marking is located on the type plate of the machinery or the optional accessory, respectively.

Unauthorized modifications of the machinery or the optional accessory will cancel this declaration.

11 Parts and Tools

NOTICE

Risk of property damage: Spare parts and consumables

Using spare parts or consumables that are not provided by ZEISS can lead to property damage.

- ▶ Only genuine spare parts and consumables supplied by ZEISS are to be used in servicing the microscope.
- ▶ Contact your ZEISS service representative for information regarding how to order spare parts and consumables.
- ▶ Unless otherwise authorized by ZEISS, all spare parts and consumables must be installed by a ZEISS service representative.

11.1 Consumables

| Required Parts/Tools | Part Number |
|---|-----------------|
| Schottky field emitter (gun) by DENKA | 000000-0302-460 |
| Aperture (single hole), 35 µm thin film | 348520-0229-000 |
| Anode aperture, 55 µm (40 nA High Resolution configuration) | 348520-0621-000 |
| Anode aperture, 90 µm (100 nA High Current configuration) | 348520-0580-001 |
| Extractor aperture | 348520-0097-001 |
| Anode aluminum seal | 348520-0266-001 |
| Copper seal gun head (single use) | 000000-0546-290 |
| Tip seal for pre-vacuum pump BOC Edwards XDS 10 | 000000-0113-989 |
| Scintillator for SE detector | 348306-8142-000 |
| Ion source (gallium) | 360100-0000-540 |
| With GIS upgrade: Precursors | On request |

11.2 Spare Parts

| Required Parts/Tools | Part Number |
|----------------------|-----------------|
| Chamber door O-ring | 000000-0476-960 |

11.3 Tools and Accessories

| Required Parts/Tools | Part Number |
|----------------------|-----------------|
| Faraday cup | 348342-8055-000 |
| 3 mm Allen key | 000000-0015-247 |

| Required Parts/Tools | Part Number |
|-----------------------------|-----------------------------------|
| 1.5 mm Allen key | 000000-0151-883 |
| Small pliers | – |
| Specimen holders | Refer to specimen holder catalog. |
| Stubs | – |
| Tweezers | – |
| Lint-free cloth | – |
| Lint-free gloves | – |

Index

A

| | |
|-------------------|-------------|
| aBSD detector | 47, 86 |
| Accessory | 131 |
| Adjusting | |
| Probe current | 80 |
| Analytic gun mode | 79 |
| Aperture | 94 |
| aSTEM detector | 50, 88, 105 |

B

| | |
|-------------------------|--------|
| Backscattered electron | 37 |
| Bakeout | 102 |
| Baking out the gun head | 102 |
| Beam mode | 33 |
| BSD detector | 46, 86 |
| BSD4 detector | 105 |

C

| | |
|---------------------------------|--------|
| CAN communication | 99 |
| Cathodoluminescence | 37, 52 |
| CCD camera | 40 |
| Chamber CCD camera | 40 |
| Chamber door seal | 101 |
| Chamberscope | 40 |
| Checking | |
| CAN communication | 99 |
| Lifetime of the ion source | 106 |
| Position of the circuit breaker | 104 |
| Temperature | 102 |
| Chemicals | 94 |
| CL detector | 52, 91 |
| Column mode | 82 |
| Consumable | 130 |
| Consumable and chemical | 94 |
| Consumables | 94 |
| Control panel | 53 |
| Cover panel, protective | 18 |

D

| | |
|---------------------------|-------------|
| Declaration of conformity | 129 |
| Detector | 38, 82 |
| aBSD | 47, 86 |
| aSTEM | 50, 88, 105 |
| BSD | 46, 86 |
| BSD4 | 105 |
| CL | 52, 91 |
| EsB | 45, 84 |
| Everhart-Thornley | 40 |
| InLens SE | 39, 83 |
| SE | 40, 83 |
| SESI | 44, 85 |
| STEM | 88 |

| | |
|---------------------|-----|
| Detector mechanics | |
| Lubricating the rod | 105 |
| Detector settings | 83 |
| Disposal | 108 |
| Microscope | 110 |
| Solid waste | 110 |
| Dual joystick | 43 |

E

| | |
|----------------------------|--------|
| Electron optical column | 31 |
| Emergency off (EMO) option | 19 |
| Emergency shutdown | 93 |
| Emitter life | 94 |
| EsB detector | 45, 84 |
| Everhart-Thornley detector | 40 |
| Extractor voltage | 81 |

F

| | |
|--------------|-----|
| Faraday cup | 80 |
| FIB column | 35 |
| Focus wobble | 76 |
| Fuse | 104 |

G

| | |
|-------------------|-----|
| Gemini II column | 31 |
| GUI | 55 |
| Gun head, bakeout | 102 |
| Gun mode | |
| Analytic | 79 |
| Imaging | 79 |
| Normal | 79 |

H

| | |
|-----------------------|----|
| Hazard | 15 |
| Biological hazard | 15 |
| Burn hazard | 15 |
| Chemical hazard | 15 |
| Electrical hazard | 16 |
| High pressure hazard | 16 |
| Magnetic field hazard | 16 |
| Mechanical hazard | 16 |
| Radiation hazard | 17 |
| Suffocation hazard | 17 |

I

| | |
|------------------------|--------|
| Image | |
| Acquisition | 68, 74 |
| Optimization | 76 |
| Saving | 78 |
| Imaging | |
| Imaging mode | 36 |
| Imaging gun mode | 79 |
| Initializing the stage | 99 |

| | | | |
|--|--------|---|---------|
| InLens SE detector | 39, 83 | Probe current | 80 |
| Installation | 64 | Calibration | 103 |
| Intended use | 13 | | |
| Ion source | | Q | |
| Checking the lifetime | 106 | Quiet Mode | 30 |
| Regenerating | 106 | | |
| Ion-sculptor focused ion beam (FIB) column | 35 | R | |
| | | Regenerating the ion source | 106 |
| J | | Related documents | 11 |
| Joystick TV angle | 100 | Replacing the chamber door seal | 101 |
| | | Requirements for personnel | 14 |
| K | | Resetting the touch alarm | 101 |
| Keyboard Shortcut | 66 | | |
| | | S | |
| L | | Safe operation condition | 14 |
| License | | Safety | 13 |
| REDUCED | 77 | Safety equipment | 18 |
| Locking device | 20 | Safety label | 22 |
| Lockout/tagout | 93 | Safety procedure | 14 |
| Lockout/tagout equipment | 21 | Saving images | 78 |
| LOTO | 21, 93 | Scintillator | 94 |
| | | SE detector | 40, 83 |
| M | | Secondary electron | 37 |
| Main Switch | 19 | Secondary electrons secondary ions detector | 44 |
| Maintenance | | Service | 14 |
| Consumable and chemical | 94 | SESI detector | 44, 85 |
| Interval | 94 | Setting the probe current | 80 |
| Procedure | 94 | Shortcuts | 66 |
| Microscope | | Shutdown | |
| Disposal | 110 | Emergency | 93 |
| Shutdown | 108 | Finishing the work session | 91 |
| Switching on | 65, 66 | Microscope | 108 |
| | | Signal detection | 37 |
| N | | SmartFIB Program Suite | 63 |
| Normal gun mode | 79 | SmartFIB User Interface | 61 |
| | | SmartSEM | |
| O | | User interface | 55 |
| Online help | 66 | SmartSEM program suite | 59 |
| Operator | | SmartSEM user interface | 55 |
| Advanced | 14 | Software | 55 |
| Basic | 14 | Software manual | 11 |
| Training | 14 | Solid waste, disposal | 110 |
| Optimizing the image | 76 | Spare part | 15, 130 |
| | | Special Keys | 66 |
| P | | Specimen current | 80 |
| Prevention of hazards | 15 | Stage | |
| Preventive maintenance | 14 | Post initialization position | 100 |
| Primary electron | 37 | Stage control | 43 |
| | | Stage initialization | 99 |
| | | Starting SmartSEM | 66 |
| | | Starting the microscope | 66 |
| | | STEM detector | 88 |
| | | Supervisor privilege | 106 |
| | | Switching off | |
| | | Gun | 102 |

| | |
|--------------|--------|
| Switching on | |
| EHT | 74 |
| Gun | 73 |
| Microscope | 65, 66 |

T

| | |
|----------------------|--------|
| Temperature | 102 |
| Tool | 131 |
| Touch alarm, reset | 101 |
| Transmitted electron | 37 |
| Troubleshooting | 96 |
| Type plate | 25, 26 |

U

| | |
|-------------------|----|
| User access level | 57 |
| User interface | 55 |
| User privilege | 57 |

V

| | |
|---------------|----|
| Vacuum system | 29 |
|---------------|----|

Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena
Germany

phone: +49 1803 33 63 34
fax: +49 3641 64 3439

info.microscopy.de@zeiss.com
www.zeiss.com/microscopy